

**DGH 6000 (SCANMATE-A)  
ULTRASONIC A-SCAN**



**OPERATOR'S MANUAL**

**For Use with Scanmate Software v4.1.x**

Equipment Manufactured By

**DGH TECHNOLOGY, INC.**



110 SUMMIT DRIVE  
SUITE B  
EXTON, PA 19341  
USA (610) 594-9100

Authorized Representative

**EMERGO EUROPE**



Prinsessegracht 20  
2514 AP, The Hague  
The Netherlands

**CE 0120**

This page intentionally left blank.

## TABLE OF CONTENTS

1	GENERAL DEVICE DESCRIPTION.....	6
2	DEVICE CLASSIFICATION.....	7
3	INDICATIONS FOR USE.....	7
4	WARNINGS AND CAUTIONS .....	7
	4.1 Meaning of Signal Words.....	7
	4.2 Device Labels .....	8
	4.3 Description of Symbols .....	8
	4.4 General Cautions and Warnings .....	9
5	PRESCRIPTION DEVICE STATEMENT .....	10
6	OPERATOR QUALIFICATIONS .....	10
7	USE OF ULTRASOUND IN OPHTHALMIC MEASUREMENTS.....	10
	7.1 Introduction to Ultrasound.....	10
	7.2 Using Ultrasonic Signals to Ascertain Correct Probe Alignment.....	11
	7.3 Ultrasonic Measurement.....	13
8	ULTRASONIC EXPOSURE AND INTENSITIES .....	15
	8.1 Tissue Exposure to Ultrasound Energy .....	15
	8.2 Ultrasonic Intensities .....	15
9	BIOMETRIC MEASUREMENT CAPABILITIES.....	16
10	INSTALLATION AND CONFIGURATION OF SCANMATE SOFTWARE.....	16
11	STARTING THE SCANMATE SOFTWARE.....	16
	11.1 Launching the Application.....	16
	11.2 Splash Screen.....	16
	11.3 Login Screen.....	17
	11.4 No USB Devices Detected Warning.....	17
	11.5 Touch Controls .....	17
12	CONFIGURING THE SCANMATE SOFTWARE.....	18
	12.1 System Preferences .....	18
	12.2 Operator Preferences.....	20
	12.3 Doctor Preferences.....	22
	12.4 Configuring Custom Velocities .....	25
	12.5 Adding IOLs to the Database.....	26
	12.6 Personalizing Lens Constants .....	28
	12.6.1 Overview.....	29
	12.6.2 Procedure .....	30
13	THE PATIENT DATA SCREEN .....	33
	13.1 Patient Data Screen Controls .....	33
	13.2 Entering a New Patient .....	34
	13.3 Searching for a Patient.....	35
	13.4 Editing a Patient Record .....	35
	13.5 Creating a New Procedure .....	36
	13.6 Entering Pre-Operative Data.....	36
	13.7 Entering Post-Operative Data .....	36
14	THE A-SCAN SCREEN.....	38
	14.1 Selecting the Operator .....	38
	14.2 Selecting the Lens and Vitreous Type .....	38

14.3	Selecting OD or OS .....	39
14.4	Measurement Modes.....	39
14.5	Automatic Measurement Mode .....	39
14.6	Manual Measurement Mode .....	40
14.7	Selecting Immersion or Contact Measurement Modes.....	42
14.8	Setting the Corneal Compression Lockout Sensitivity.....	42
14.9	Performing a Contact Measurement .....	43
14.10	Performing an Immersion Measurement .....	47
14.11	Saving a Measurement File.....	51
14.12	Saving a Video Buffer .....	52
14.13	Reviewing Measurements and Videos.....	52
14.14	Voice Controls .....	52
15	THE B-SCAN SCREEN .....	54
16	THE IOL CALCULATOR SCREEN .....	55
16.1	IOL Calculator Screen Controls .....	56
16.2	Patient Information .....	56
16.3	Selecting the Doctor.....	56
16.4	Eye Selection .....	56
16.5	Post Refractive Selection .....	57
16.6	Measurement Source.....	57
16.7	Corneal Refractive Index .....	57
16.8	Corneal Power (K1 and K2) .....	58
16.9	Desired Refraction .....	58
16.10	Selecting IOLs .....	58
16.11	Selecting IOL Power Calculation Formulas .....	58
16.12	Selecting Personalized Lens Constants.....	59
16.13	K Correction Method Selection .....	59
16.14	IOL Calculator Outputs .....	59
16.15	Saving Pre-Operative Data .....	59
16.16	Post Refractive Calculation Methods .....	60
16.17	Using the Post Refractive Calculator.....	63
16.18	Performing IOL Calculations.....	64
16.19	Overview of IOL Power Calculation Formulas .....	66
16.20	Factors Impacting the Accuracy of IOL Power Calculations .....	67
16.20.1	Axial Length .....	67
16.20.2	Corneal Power.....	67
16.20.3	Postoperative Anterior Chamber Depth.....	68
16.20.4	Surgical Technique .....	68
16.20.5	Implant Power.....	68
17	CREATING REPORTS .....	69
17.1	IOL Calculator Report .....	69
17.2	A-Scan Short Report.....	69
17.3	A-Scan Custom Report .....	69
17.4	Using Reports .....	70
17.5	Opening Reports .....	70
17.6	Exporting A-Scan PDF Reports.....	70
18	DATABASE MANAGEMENT.....	71

18.1	Importing and Exporting Data .....	71
19	ELECTROMAGNETIC COMPATIBILITY.....	72
20	CARE AND MAINTENANCE .....	73
20.1	Care of Transducer.....	73
20.2	Cleaning and Disinfecting the Transducer Tip.....	74
20.3	Cleaning and Disinfecting the Immersion Shell.....	77
20.4	Care of USB Interface Module .....	77
20.5	Cleaning the USB Interface Module.....	77
20.6	Operating Conditions.....	77
20.7	Verifying A-Scan Calibration.....	77
20.8	Storage .....	78
20.9	Transportation.....	79
20.10	Disposal .....	79
21	TROUBLESHOOTING.....	80
21.1	Unable to Get a Measurement .....	80
22	WARRANTY.....	80
23	LIFETIME / SHELF-LIFE.....	81
24	CUSTOMER SERVICE .....	81
APPENDIX A	COMPUTER SYSTEM SPECIFICATIONS .....	82
APPENDIX B	SCANMATE-A SPECIFICATIONS .....	83
APPENDIX C	REFERENCES .....	84
APPENDIX D	TERMS AND ABBREVIATIONS .....	86

## **1 General Device Description**

The DGH 6000 is a diagnostic ultrasound device used by professionals in the ophthalmic field to perform axial length (AXL), anterior chamber depth (ACD), and lens thickness (LT) measurements of the human eye. Measurements are obtained using ultrasonic pulse echo technology, whereby short bursts of ultrasonic energy are transmitted and the resulting echoes are captured, amplified, filtered and processed.

The DGH 6000 is comprised of three main components: a detachable 10-MHz transducer, a USB interface module and the Scanmate software application.

The transducer uses a single element piezoceramic with a starting frequency of 12.5 MHz +/- 1.25 MHz. This ceramic, after it has been damped, runs at a center frequency of 10.0 MHz nominal. The tip of the transducer has a spherical radius of curvature, which focuses the acoustic beam at 23.0 mm nominal. It is also equipped with a fixation LED to facilitate probe alignment.

The USB Interface Module contains the electronics used to pulse the transducer as well as to measure, filter and amplify the resulting echoes from internal eye structures. The USB Interface Module is powered and controlled via a USB 2.0 cable, which must be connected to a PC or laptop computer running the Scanmate software application.

The Scanmate software interprets the digital signals transmitted by the USB Interface Module and displays a live “Amplitude Scan” which shows the relative magnitude and distance of the echo spikes received by the transducer. Specific time distances between the captured echo spikes are then measured and converted into distance information. A digital signal processing algorithm is used to analyze the waveform, ensuring that the probe is properly aligned at the time of measurement. The software assists the user in performing replacement IOL power calculations using established formulas. The user can save videos or measurement files or create reports to document the results of the exam.

The captured A-Scan waveforms, associated measurements, and IOL calculations should be evaluated only by trained medical personnel. The DGH 6000 does not conduct nor provide any diagnoses, drive therapeutic outputs, or determine or maintain any life sustaining or life threatening equipment.

## 2 Device Classification

Device: System, Imaging, Pulsed Echo, Ultrasonic  
Panel: Radiology  
Product Code: IYO  
Device Class: II  
Regulation Number: 21 CFR 892.1560

Device: Diagnostic Ultrasonic Transducer  
Panel: Radiology  
Product Code: ITX  
Device Class: II  
Regulation Number: 21 CFR 892.1570

## 3 Indications for Use

The intended use of the DGH 6000 is the measurement of axial length, anterior chamber depth, and lens thickness of the human eye. The DGH 6000 is also intended to calculate the optical power of an IOL that is to be implanted during cataract surgery. The DGH 6000 is intended to be used solely by qualified medical professionals. Clinical consideration and judgment should be applied when using the DGH 6000.

## 4 Warnings and Cautions


### 4.1 Meaning of Signal Words

In this manual, the signal words “Warning” and “Caution” are used to highlight important safety and usage instructions. All users of the DGH 6000 must understand the meanings of these signal words.

Signal Word	Meaning
 <b>WARNING</b>	Indicates a potentially hazardous situation which if not avoided could cause injury or harm to the equipment or erroneous results.
 <b>CAUTION</b>	Indicates a potentially hazardous situation which if not avoided may result in minor injury or harm to the equipment.

## 4.2 Device Labels

The label shown below is located on the underside of the USB Interface Module.



## 4.3 Description of Symbols



This symbol indicates the degree of protection against electric shock. The DGH 6000 Scanmate-A is classified as type BF equipment.



This symbol indicates that the DGH 6000 Scanmate-A is an IEC Protection Class II (double insulated) device.



This symbol instructs the operator to read the operating manual.



This mark indicates that Notified Body 0120 (SGS United Kingdom Ltd) has certified the management system of DGH Technology, Inc. meets the requirements applicable requirements of 21 CFR 1010 (Performance Standards for Electronic Products: General) and 21 CFR 1050 (Performance Standards for Sonic, Infrasonic, and Ultrasonic Radiation-Emitting Products). The device also conforms to the following International Standards:

- EN 60601-1: Medical electrical equipment – Part 1: General requirements for safety – IEC 60601-1
- EN 60601-1-2: Medical electrical equipment – Part 1: General requirements for safety. Collateral standard: Electromagnetic compatibility requirements and tests. IEC 60601-1-2
- NEMA Standard Publication UD-2: Acoustic Output Measurement Standard for Diagnostic Ultrasound Equipment
- NEMA Standard Publication UD-3: Standard for Real-Time Display of Thermal and Mechanical Acoustic Output Indices on Diagnostic Ultrasound Equipment



This symbol indicates that Emergo Europe is the European Authorized Representative for this device.



This symbol indicates that DGH Technology, Inc. is the manufacturer of the DGH 6000 Scanmate-A device. The YYYY under the symbol indicates the year the device was manufactured.

**REF** This symbol indicates that the model number of this device is DGH 6000.

**SN** This symbol indicates the serial number of the device. YYYY indicates the year of manufacture and XXXX indicates the unit number.



This symbol located on the DGH 6000 indicates that the equipment consists of electronic assemblies and other components that may be subject to Directives 2002/96/EC, 2003/108/EC, and 2002/95/EC of the European parliament, which advises that electrical and electronic devices must not be disposed of as normal domestic refuse. In order to prevent environmental risks or endangerments by non-professional disposal, the disposal of this product, including any accessories, must comply with valid practices as outlined in Directives 2002/96/EC, 2003/108/EC, and 2002/95/EC and local regulations. All electronic components and systems should be returned to Original Manufacturer for disposal.

#### 4.4 General Cautions and Warnings



#### CAUTION

Transducers must be cleaned after each use. Cleaning the transducer is an essential step prior to effective disinfection. Follow the manufacturer’s instructions when using disinfectants.



#### WARNING

Do not allow sharp objects, such as scalpels or cauterizing knives, to touch the transducer or cables.



#### WARNING

Equipment not suitable for use in the presence of flammable mixtures.

## 5 Prescription Device Statement



The DGH 6000 is a prescription device and is only to be used by, or under the supervision of, a licensed physician.

## 6 Operator Qualifications

This DGH 6000 is intended to be used by trained medical professionals. The medical professional operating the DGH 6000 must have a general knowledge of the use of ultrasonic medical devices. Use of the DGH 6000 requires adequate dexterity to position the probe safely. The DGH 6000 uses audio feedback to inform the operator of the scan status.

## 7 Use of Ultrasound in Ophthalmic Measurements

### 7.1 Introduction to Ultrasound

Ultrasound offers a non-invasive method to examine the interior of solid objects. Ultrasonic pulses consist of sound waves of a frequency level too high to be heard by the human ear. When a sound impulse strikes an interface, some sound is reflected, and some sound is transmitted. Because some sound will pass through the surface and be reflected by the next surface, complex structures can be examined with ultrasound. When ultrasound penetrates an object with several interfaces, the reflected ultrasound can be observed on a graphic display as a two-dimensional waveform with peaks that are related to the positions of the interfaces.

The DGH 6000 transducer emits ultrasound pulses and detects ultrasound signals that have been reflected back. The time delay between the echoes is used to calculate distances between surfaces in the eye.

---


**Note:** Ultrasound cannot travel through air because air is not dense enough for the high frequency waves to propagate. Ultrasonic measurements must therefore be performed by direct contact or through a denser medium such as water.

---


## 7.2 Using Ultrasonic Signals to Ascertain Correct Probe Alignment

Sound travels in straight lines, so the direction of reflected sound is based solely on its angle of incidence. Sound hitting an interface perpendicularly will reflect back along the same path that it approached (Figure 7-A). Sound hitting an interface at an angle will reflect at an angle away from the source (Figure 7-B). The transmitted sound will continue on at a lesser amplitude because of reflected energy lost at the interface.

When reflected ultrasound is shown as a two-dimensional waveform, the peaks are related to the positions of the interfaces. By comparing the relative height (intensity) of the peaks, one can determine the angle at which the sound is striking it (Figures I-A and I-B). Steadily diminishing peaks are an indicator that the ultrasound is not perpendicular to the interfaces.




**Figure 7-A** *Sound hitting an interface perpendicularly.*




**Figure 7-B** *Sound hitting an interface at an angle.*

Using these properties of ultrasound, the alignment of an ultrasound beam through the eye can be determined. Proper alignment is crucial to the accuracy of measurements that will be used for IOL calculations. Figure 7-C illustrates an ultrasonic pattern typical when alignment along the visual axis is met. Please note the two high, clear peaks representing the anterior lens and posterior lens interfaces, along with a strong peak representing the retinal interface. Figure 7-D illustrates a pattern representing misalignment. Incorrect alignment will result in improper measurements.

**Figure 7-C** Correct alignment along the visual axis.




**Figure 7-D** Incorrect alignment along the visual axis.



The DGH 6000 Scanmate-A incorporates a pattern recognition program that automatically checks for proper alignment. It first looks for the proper pattern of reflectance peaks from each expected interface. Then the software examines the retinal spike for particular signal characteristics which are produced only by retinal interfaces. Incorrect signal patterns will not be measured until the alignment is corrected. Audible feedback will guide the operator to proper alignment.

### 7.3 Ultrasonic Measurement

The speed of sound increases in denser materials. Liquids or substances containing large amounts of water conduct ultrasound very well; air does not conduct ultrasound. Figure 7-E shows a list of the established speeds of sound for ophthalmic structures.



**Figure 7-E** *Speeds of sound for ophthalmic structures.*

Using the relationship between the density of a material and the speed of sound, ophthalmic A-Scans obtain distances in the eye by performing a two-step process. First, a pulse of sound is timed as it travels through the eye, reflects off the retina, and returns to the transducer. Second, the lengths are calculated based on the travel time and the speed of sound through the eye:

$$\text{distance} = \frac{\text{velocity} \times \text{time}}{2}$$

Early ophthalmic ultrasonic units used an average speed of sound of 1550 m/s for phakic eyes and 1532 m/s for aphakic eyes. The DGH 6000 Scanmate-A lets the user select the type of eye being examined, and will then automatically use the proper speed for each part of the eye. Using specific velocities for each part of the eye increases the accuracy of measurements. Ex:

$$\text{lens thickness (normal)} = \frac{1641 \text{ m/s} \times \text{time (s)}}{2}$$

Cornea thickness in all eye types is assumed to be 550 μm, with a velocity of 1641 m/s. This assumption is part of the ACD calculation. Anterior Chamber Depth is defined as the measured distance from anterior vertex of the cornea and the anterior vertex of the lens.

In pseudophakic eyes, the lens thickness is not measured. A standard thickness is assumed, depending on the type of implant present. The following table summarizes the assumed thicknesses and velocities.

Material	Velocity (m/s)	Thickness (µm)
PMMA IOL	2780	700
Silicone IOL	980	1350
Acrylic IOL	2180	800

**Note:** The velocities table within the Scanmate application displays the preset velocities and allows the user to specify a custom velocity and thickness to be used for calculation. Custom velocities are selected as Custom Aphakic, Custom Phakic, or Custom Pseudophakic lens types.

**Standard Velocities**

	Speed (m/s)	Thickness (mm)
Cornea	1641.00	Measured
Aqueous	1532.00	Measured
<b>Lens</b>		
Cataract	1629.00	Measured
Dense Cataract	1590.00	Measured
Normal	1641.00	Measured
PMMA IOL	2780.00	0.70
Silicone IOL	980.00	1.35
Acrylic IOL	2180.00	0.80
<b>Vitreous</b>		
Normal	1532.00	Measured
Silicone Oil 1000cs	980.00	Measured
Silicone Oil 5000cs	1040.00	Measured

Aphakic  
 Phakic  
 Pseudophakic

Custom Velocities	Speed (m/s)	Thickness (mm)
Cornea	1641.00	0.550
Aqueous	N / A	N / A
Lens	N / A	N / A
Aphakic	1532.00	Measured

NOTE: Changes to custom velocities do not take effect until the lens type is re-selected. This can be caused through explicit selection, starting a new measurement, or loading a video.

## 8 Ultrasonic Exposure and Intensities

### 8.1 Tissue Exposure to Ultrasound Energy

The ultrasound energy emitted by the DGH 6000 is low intensity and will have no adverse effects on the patient and/or operator. However, the operator is still cautioned to perform examinations using the principle of **ALARA** (As Low As Reasonably Achievable). All examinations should be done so that the patient receives as little ultrasound radiation as possible. Do not hold the probe against the eye or other tissue with the system activated except when making a measurement. Do not make unnecessary measurements.

### 8.2 Ultrasonic Intensities

The DGH 6000 has only one mode, and ultrasonic intensity settings are not under the control of the operator. Thus, the values below are the values to be expected for a typical transducer.

Since the DGH 6000 is not capable of exceeding either a TI of 1.0 or an MI of 1.0 in any operating mode, the output of the system is reported as shown in the Table below.

The appropriate Thermal Index is the Thermal Index for Soft Tissue, TIS, for the non-scanning case with a beam aperture of less than 1.0 cm.

#### Output Summary Table

Transducer Model (used with DGH 6000)	$I_{\text{spta.3}}$	TI Type	TI Value	MI	$I_{\text{pa.3}} @ \text{MI}_{\text{max}}$
DGH6006	0.469 mW/cm <sup>2</sup>	TIS non-scan, $A_{\text{aprt}} < 1.0$	0.0004	0.131	3.49 W/cm <sup>2</sup>

The acoustic output values given above are based on a presumed attenuation of ultrasound on tissue, as developed by the U.S. Food and Drug Administration in 1985, and later incorporated into other international standards.

The attenuated intensity in the eye at the transducer focus (corresponding to maximum intensity) may be calculated according to the formula recommended by the FDA:

$$I_t = I_w \times e^{(-0.069 \times f \times z)}$$

where  $I_t$  is the estimated in situ intensity,  $I_w$  is the measured intensity in water at the focus of the transducer,  $f$  is the ultrasonic frequency, and  $z$  is the distance from the face of the probe to the transducer focus, which is the point of measurement (23 millimeters).

The nominal piezoceramic (crystal) frequency of these transducers is 10 MHz. The actual frequency of a particular transducer may vary from this value. The tissue calculations above were done with the measured frequency of the transducer used for the tests.

## 9 Biometric Measurement Capabilities

The following table shows the measurement range and capabilities for the DGH 6000 (Scanmate A).

Parameter	Range	Accuracy	Resolution
Axial length	15 mm to 40 mm	+/- 0.1mm	0.01 mm
Anterior chamber depth	2.0 mm to 6.0 mm	+/- 0.1 mm	0.01 mm
Lens thickness	2.0 mm to 7.5 mm	+/- 0.1 mm	0.01 mm

## 10 Installation and Configuration of Scanmate Software

Refer to the Scanmate Installation Manual for information on installing and configuring the software.

## 11 Starting the Scanmate Software


### 11.1 Launching the Application



Once installed, the “Scanmate” shortcut appears on the Windows desktop and in the start menu. Click on the desktop icon to start the DGH Scanmate application.

### 11.2 Splash Screen

The Scanmate splash screen will appear while the application loads.



### 11.3 Login Screen

A single username and password is used to gain access to the software and database for all users. By default, the Scanmate software is set to automatically log in when the application is started using the username and password specified in **System Preferences**. To change this setting, uncheck the “Automatic Login” preference in the **System Preferences** menu. If a login is requested, enter the username and password created during software installation.



### 11.4 No USB Devices Detected Warning

If the USB Interface Module is not attached, a warning box will appear.



Clicking “OK” will complete the log-in and allow the Scanmate software to be used without the USB Interface Module. Although no scans can be completed, the software can still be used to perform IOL calculations or view completed reports, saved measurements and saved videos.

If the Scanmate software is being used without the USB Interface Module attached, the software will require probe key authentication before using, and after every 20 hours of use. Warnings will appear every hour after 15 hours of use. To provide authentication, plug a USB Interface Module into a USB port; authentication will complete in a few seconds.


### 11.5 Touch Controls

The Scanmate software can be operated using touch controls on systems that have a touch-capable display. Buttons, sliders and combo boxes can be operated by touching the screen.

## 12 Configuring the Scanmate Software

### 12.1 System Preferences

The “System Preferences” window provides controls to set various configuration items for the system. The System Preferences window can be accessed by selecting **Preferences** → **System** from the Menu Bar.



The screenshot displays the System Preferences window with two main sections: System Options and Login Configuration. The System Options section includes fields for Practice Name (DGH Technology, Inc.), Practice Address Line 1 (110 Summit DR; Suite B), Practice Address Line 2 (Exton, PA 19341), Date Format (MM/dd/yyyy), Report PDF Export Location, Default Keratometer Index (nc) (1.3375), Enable Voice Control (checked), Hey Flex, and Voice Control Sensitivity. The Login Configuration section includes Database Type (SQL Server), Database Location ((local)\DGHSERVER), Database Name (DGH-Scanmate), Username (User), Password (masked), Hint (!DGHTest1), and Automatic Login (checked). A Change Configuration button is located at the bottom right of the Login Configuration section. At the bottom of the window, there is a caution message: "Caution: Accepting changes to System Preferences may alter currently loaded data. Please review data on all forms." and OK and Cancel buttons.

The **Practice Name**, **Practice Address Line 1** and **Practice Address Line 2** textboxes allow the user to enter practice information that will appear on reports generated by the software.

The **Date Format** textbox allows the user to select either MM/dd/yyyy or dd/MM/yyyy date format for the software.

The **Automatic Login** checkbox allows the user to select if login is required upon system startup or if the user is automatically logged in. By default, this checkbox is enabled, indicating that the username and password will be entered automatically.

The **Default Keratometer Index (nc)** box is provided to give the user the ability to configure the default index of refraction used on the IOL Calculator screen.

 **WARNING**

To avoid IOL power calculation errors, it is important that the corneal refractive index (nc) of the keratometer used to calculate corneal power matches the default refractive index (nc) set in System Preferences for the DGH Scanmate.

The default nc should be set before K1 and K2 are entered for a patient. The nc associated with each set of corneal power measurements can be reviewed or adjusted on the **IOL Calculator** Screen. Failure to adjust the default nc to match the nc of the keratometer used to measure corneal power can result in errors in the IOL Calculation. See Section 16.20 of this guide for more details.

The **Report PDF Export Location** setting allows the user to specify the default directory where pdf reports will be exported. The default directory can be any local or mapped network drive.

The **Enable Voice Control** checkbox allows the user to enable or disable voice control. If voice control is enabled, the user can select how to initiate voice commands from the textbox below. Voice commands may be configured to be initiated by pressing the Enter key or by a verbal command of “Hey Flex”. Voice control sensitivity can be adjusted in this menu as well.

The **Database Type** textbox lists the database type currently being used by the Scanmate software. For more information on database types, refer to the Scanmate Installation Guide.

The **Database Location** textbox specifies the location of the DGH Database Server hosting the DGH-Scanmate database. To change the Database Location for the DGH Scanmate application, select the “Change Configuration” button and follow the prompts. For more information on configuration of the software, refer to the Scanmate Installation Guide.

The **Database Name** textbox displays the name of the DGH-Scanmate database. The default database name is created automatically when installing the system and is not user configurable.


The **Username** textbox displays the username that the Scanmate application uses to connect to the database. To change the Username for the DGH Scanmate application, select the “Change Configuration” button and follow the prompts. For more information on configuration of the software, refer to the Scanmate Installation Guide.

The **Password** textbox displays the password that the Scanmate application uses to connect to the database. To change the Hint that the DGH Scanmate application displays, select the “Change Configuration” button and follow the prompts. For more information on configuration of the software, refer to the Scanmate Installation Guide.

The **Hint** textbox allows the user to enter a hint for the password. If a user is logging in manually, clicking the “Show Hint” button at log-in will display the hint text. To change the Hint that the DGH Scanmate application displays, select the “Change Configuration” button and follow the prompts. For more information on configuration of the software, refer to the Scanmate Installation Guide.

## 12.2 Operator Preferences

The “Operator Preferences” window provides controls to identify Operators and configure their preferences. The Operator Preferences window can be accessed by selecting **Preferences** → **Operator** from the Menu Bar.



---

**Note:** Refer to the DGH 8000 Scanmate-B or the Scanmate Flex User Guide for information on configuring B-Scan related Operator Preferences.


---

The **New** button allows the user to create new operators and assign preferences to the newly created operator.

The **Edit Name** button allows the user to make changes to an existing operator’s preferences. To change an operator’s preferences, click the “Edit” button while the desired operator’s name is selected in the “Name” dropdown box. Changes to operator

preferences are saved by clicking the “Apply” button, or can be discarded by clicking “Cancel”. The “OK” button will save the current settings and close the “Operator Preferences” window.


The **Delete** button will delete an operator from the system. When clicked, the user is prompted for confirmation prior to deleting the operator’s preferences.



The **Name** combo box selects an existing operator from the database to view, edit, or delete their preferences. Selecting an operator’s name will display that operator’s preferences on the screen. The **Default Operator** checkbox selects which operator’s preferences will be automatically loaded when the software is started.



The **Notes** textbox allows the user to associate notes to an operator’s preferences.



The A-Scan **Video Playback Speed** slide bar adjusts the frames per second (fps) that the A-Scan video is played for the currently selected operator. The speed of video playback can be adjusted from 1 to 20 fps.

The A-Scan **Video Buffer Frames** slide bar adjusts the size of the buffer used when capturing A-Scan video. Moving the slide bar from left to right increases the number of frames stored in the buffer. Increasing the number of frames that can be stored in the buffer allows the user to determine the amount of video that is captured before recording is stopped. The video buffer can be configured to store from 0 (zero) to 1024 frames.

The A-Scan **Probe LED Intensity** slide bar adjusts the default intensity of the LED on the probe for the currently selected operator. The intensity (brightness) can be adjusted from 0 (zero) to 100.

### 12.3 Doctor Preferences

The “Doctor Preferences” window provides controls for the user to identify Doctors and configure their preferences. Doctor Preferences define the default protocol for performing A-Scans and IOL Calculations. The Doctor Preferences window can be accessed by selecting **Preferences → Doctor** from the Menu Bar.

The screenshot shows the 'Doctor Preferences' window with the following details:

- Doctor Section:** Includes 'New', 'Edit Name', and 'Delete' buttons. A 'Name' dropdown menu is set to 'Default Preferences'. A 'Default Doctor' checkbox is checked. A 'Notes' text area contains 'Default Preferences'.
- A-Scan Auto Stop Scan Criteria:** Three sliders for 'AXL Std. Deviation Threshold', 'ACD Std. Deviation Threshold', and 'LT Std. Deviation Threshold', all set to 0.05 mm.
- A-Scan Default Measurement Mode:** Radio buttons for 'Contact' and 'Immersion', with 'Immersion' selected.
- A-Scan Timing:** Sliders for 'Measurement Delay' (0.0 sec) and 'Alignment Timeout' (20.0 sec).
- IOL Calculator Default Settings:** 'Min Cases' slider set to 25. 'Meas Source' dropdown set to 'Unit Average'. Three tabs for 'IOL 1', 'IOL 2', and 'IOL 3'. Below the tabs are dropdowns for 'Mfg:', 'Model:', and 'Formula:'. Buttons for 'A-Scan Velocities', 'IOL Configurations', and 'Personalized Lens Constants' are at the bottom.
- Footer:** A caution message: 'Caution: Accepting changes to Doctor Preferences may alter currently loaded data. Please review data on all forms.' and 'OK', 'Cancel', and 'Apply' buttons.

The **New** button allows the user to create new doctor profiles and assign preferences to the newly created doctor.

The **Edit Name** button allows the user to make changes to an existing doctor’s preferences. The user can change a doctor’s preferences by clicking this button while the desired doctor’s name is selected in the “Name” combo box.


The **Delete** button deletes a doctor from the system. When clicked, the user is prompted for confirmation prior to deleting the doctor’s preferences.

This close-up shows the 'Name' dropdown menu with 'Default Preferences' selected and the 'Default Doctor' checkbox checked.

The **Name** combo box selects an existing doctor from the database to view, edit, or delete their preferences. The **Default Doctor** checkbox selects which doctor will be automatically placed on the Patient Data and IOL Calculator Screens. If only one doctor exists, that doctor is automatically set as the default doctor.

This close-up shows the 'Notes' text area containing the text 'Default Preferences'.


The **Notes** text box is provided to allow the user to associate notes to a doctor’s preferences.




The **AXL Std. Deviation Threshold** textbox allows the user to enter the AXL standard deviation threshold for the currently selected doctor. This threshold defines the standard deviation that the eight (8) axial length measurements in the measurement bank must fall within for the system to automatically stop scanning. See Sections 14.4, 14.5 and 14.6 for more details on automatic and manual measurement modes.

The **ACD Std. Deviation Threshold** textbox allows the user to enter the ACD standard deviation threshold for the currently selected doctor. This threshold defines the standard deviation that the eight (8) anterior chamber depth measurements in the measurement bank must fall within for the system to automatically stop scanning. See Sections 14.4, 14.5 and 14.6 for more details on automatic and manual measurement modes.

The **LT Std. Deviation Threshold** textbox allows the user to enter the lens standard deviation threshold for the currently selected doctor. This threshold defines the standard deviation that the eight (8) lens thickness measurements in the measurement bank must meet for the system to automatically stop scanning. See Sections 14.4, 14.5 and 14.6 for more details on automatic and manual measurement modes.



The **A-Scan Default Measurement Mode** radio buttons allow the user to select either contact or immersion as the default A-Scan measurement mode for the currently selected doctor.



Variables in timing used while aligning the probe can be adjusted in the **A-Scan Timing** box. Setting **Measurement Delay** to 0 ms will make the DGH 6000 take measurements as soon as alignment is achieved. Setting some Measurement Delay is useful in ensuring that a variety of different measurements are used to create an average.

The **Alignment Timeout** setting allows the user to specify how long the A-Scan probe can be appanated before a message appears indicating why the waveforms aren't meeting the measurement criteria.

Min Cases:  Before Automatically Updating Personalized Lens Constants (When Automatically Update Constants is Selected).

Meas Source:

The **Min Cases Before Using Personalized Constants** textbox allows the user to set the minimum number of Procedures that must be performed before a Personalized Lens Constant is calculated for a given IOL model. See Section 12.6 for more details on calculating personalized lens constants.

The **Meas Source** combo box allows the user to choose the measurement source that the IOL Calculator uses by default for the currently selected doctor. The user can choose to use DGH 6000 Average, DGH 6000 Single, or User Input as the selected doctor's preferred default measurement source. The IOL Calculator Screen will automatically load using the Default Measurement Source for the selected doctor. See Section 16.6 for more information on selecting a Measurement Source.

IOL 1 | IOL 2 | IOL 3

Mfg:

Model:

Formula:

The **IOL 1 / IOL 2 / IOL 3** tabbed display allows the user to specify up to three preferred IOLs for the currently selected doctor. The default IOL manufacturer (Mfg), model, and formula are selected from searchable combo boxes for each preferred IOL. The IOL Calculator Screen will automatically load the Preferred IOL tabs for the currently selected doctor. IOLs are created in the IOL Information Table (see Section 12.5). If more than three IOLs exist for a doctor, the additional IOLs can be selected on the IOL Calculator Screen.

## 12.4 Configuring Custom Velocities

### A-Scan Velocities

The **A-Scan Velocities** button on the “Doctor Preferences” screen provides the user access to the Velocities Table dialog. From this screen, the user can view the speeds of sound for each lens and vitreous type. These values are used by the software in calculations interpreting measurements taken by the DGH 6000. This screen also allows the user to enter desired speeds of sound for custom phakic, aphakic, and pseudophakic lenses. See Section 7.3 for more information about using the speed of sound in ultrasonic measurement.

Standard Velocities		
	Speed (m/s)	Thickness (mm)
Cornea	1641.00	Measured
Aqueous	1532.00	Measured
<b>Lens</b>		
Cataract	1629.00	Measured
Dense Cataract	1590.00	Measured
Normal	1641.00	Measured
PMMA IOL	2780.00	0.70
Silicone IOL	980.00	1.35
Acrylic IOL	2180.00	0.80
<b>Vitreous</b>		
Normal	1532.00	Measured
Silicone Oil 1000cs	980.00	Measured
Silicone Oil 5000cs	1040.00	Measured

Custom Velocities		
	Speed (m/s)	Thickness (mm)
Cornea	1641.00	0.550
Aqueous	N / A	N / A
Lens	N / A	N / A
Aphakic	1532.00	Measured

NOTE: Changes to custom velocities do not take effect until the lens type is re-selected. This can be caused through explicit selection, starting a new measurement, or loading a video.

**Note:** The Scanmate software can only store one set of custom velocities; changing these values under one doctor’s preferences will change them for all users.

## 12.5 Adding IOLs to the Database

### IOL Configurations

The **IOL Configurations** button opens the IOL Information Table dialog. From this screen, the user can create, edit, view, or delete IOLs from the database. The IOL Information Table can only be opened when a doctor's preferences are being created or edited. IOLs are specific to the doctor that creates them. If more than one doctor will be using the DGH 6000, they will each need to enter their preferred IOLs.

The **Doctor** combo box selects an existing doctor from the database to view, edit, add, or delete IOL information that is associated with the specified doctor.

The **New** button creates new IOL records. When this button is clicked, previously read-only textboxes become enabled and editable.

The **Edit** button allows the user to make changes to an existing IOL configuration for the selected Doctor.

The **Delete** button deletes an IOL configuration for the selected doctor. The IOL information record will not be removed from the database, but it is marked as inactive.

---

**Note:** When the lens constants specified by the IOL manufacturer are changed, the personalized lens constants will be recalculated for the IOL record if “Automatic Personalized Lens Constants Update” is enabled.

---

### WARNING

Editing the manufacturer's A-constant will impact the Personalized Lens Constant Calculations for the selected IOL.

**IOL**

New Edit Delete

Manufacturer

Model

A-Constant

Step Range Type

The **Manufacturer** combo box allows the user to view, enter and edit the manufacturer of an IOL.

The **Model** combo box allows the user to view, enter and edit the model of the IOL.

The **A-Constant** is entered by the user based on the IOL manufacturer's specifications.

The **Step Range** combo box has two possible selections: ½ Diopter and ¼ Diopter.

The **Type** combo box has two possible selections: Posterior and Anterior.

**Conversions**

Calculate from Mfg. A-Const.

A-Const.	118.05	SRK-II
A-Const.	118.05	SRK-T
pACD	4.998	Hoffer Q
SF	1.25	Holladay 1
a0	1.308	Haigis
a1	0.400	Haigis
a2	0.100	Haigis
ACD Const.		Binkhorst II (Optional)

The **Calculate from Mfg. A-Const** button will become active when the Manufacturer's A-Constant has been entered on the left. Clicking this button will automatically calculate and display all constants needed for the SRK®-II, SRK®-T, Hoffer® Q, Holladay 1, and Haigis formulas. These values can also be entered individually.

The ACD Constant used for the Binkhorst II formula is not calculated from the Manufacturer's A-Constant and should be entered here if the Binkhorst II formula will be used with this IOL.

**Notes**

The **Notes** textbox allows the user to record any notes associated with an IOL configuration.

## 12.6 Personalizing Lens Constants

Personalization of lens constants helps to improve the accuracy of a given IOL power calculation formula by correcting for consistent prediction errors caused by:

- Surgical Technique
- A-Scan Biometry Measurements
- Keratometer Biometry Measurements
- IOL Style / Geometry

The **Personalized Lens Constants Table** allows the user to manage how Personalized Lens Constants are calculated. The table displays all Pre-Operative, Post-Operative, patient, and procedure information associated with the currently selected Doctor and IOL. The information in the table is placed in a default column order; however, the user is provided with the ability to rearrange the columns into any order. The ability to sort the entire table allows the user to rearrange the rows based on the selected column.

The Personalized Lens Constants Table is opened by selecting from the Menu Bar **Preferences → Doctor**, then choosing to edit the selected doctor, and then clicking the **Personalized Lens Constants** button.

Doctor:

Mfg:

Model:

A-Constant:

ACD:

**Current Personalized Constants**

A-Const.	<input type="text" value="118.05"/>	SRKT	<input type="text" value="a0"/>	<input type="text" value="1.308"/>	Haigis	Min. Cases	<input type="text" value="25"/>
pACD	<input type="text" value="4.998"/>	Hoffer Q	<input type="text" value="a1"/>	<input type="text" value="0.400"/>	Haigis	Selected Cases	<input type="text" value="0"/>
SF	<input type="text" value="1.25"/>	Holladay 1	<input type="text" value="a2"/>	<input type="text" value="0.100"/>	Haigis	<input type="button" value="Edit"/>	

Automatically Update Personalized Lens Constants (After Min. # of Cases Achieved)

**New Personalized Constants**

A-Const.	<input type="text" value="118.05"/>	SRKT	<input type="text" value="a0"/>	<input type="text" value="1.308"/>	Haigis	Selected Cases	<input type="text" value="0"/>
pACD	<input type="text" value="5.020"/>	Hoffer Q	<input type="text" value="a1"/>	<input type="text" value="0.400"/>	Haigis		
SF	<input type="text" value="1.28"/>	Holladay 1	<input type="text" value="a2"/>	<input type="text" value="0.100"/>	Haigis		

Use Record	Post Refractive Method	Patient ID	Surgery Date	Measurement Source	AXL	ACD	K1	K2	Desired Rx	Implanted Lens Power

## 12.6.1 Overview

The DGH 6000 assists the user in performing personalized lens constant calculations for the following formulas:

- SRK®-T (A-Constant)
- Hoffer® Q (pACD)
- Holladay 1 (SF)
- Haigis (Single Optimization of a0)

---

**Note:** The DGH 6000 performs a Single Optimization (a0) for the Haigis Formula. Triple Optimization (a0, a1 and a2) is unavailable at this time.

---

The Personalized Lens Constants Table contains the following features:

**Doctor** - The “Doctor” combo box selects an existing doctor from the database. Personalized Lens Constant calculations are performed for the currently selected doctor.

**Manufacturer** - The “Mfg” combo box selects the manufacturer of an IOL. The list of IOL manufacturers is limited based on the manufacturers that are associated with the selected doctor.

**Model** - The “Model” combo box selects the model of an IOL. The list of IOL models is limited based on the models that are associated with the selected doctor and the selected manufacturer.

**A-Constant** - The “A-Constant” textbox allows the user to view the A-Constant for the currently selected IOL. This is the Manufacturer’s A-Constant that was entered in the IOL Information Table.

**ACD** - The “ACD” textbox allows the user to view the ACD-Constant for the currently selected IOL. This is the ACD-Constant that was entered in the IOL Information Table.

**Current Personalized Constants** - The “Current Personalized Constants” group box displays the lens constants that are currently being used for the selected Doctor and IOL. All textboxes within this group box are read-only unless the user clicks the “Edit” button. Once the “Edit” button is selected, all textboxes become editable to allow the user to manually change any of the current constant values for the selected Doctor, IOL Manufacturer and IOL Model. If a value outside of the accepted range for each variable is entered, it will be highlighted in red and must be corrected before continuing.

**New Personalized Constants** - The “New Personalized Constants” group box displays the lens constants as they are calculated by selecting and deselecting records in the “Personalized Lens Constants” table. This allows the doctor to view the updated calculation before accepting it with the “Up Arrow” button.

**Up Arrow** – When this button is clicked, the values that are displayed within the “New Personalized Constants” textboxes are stored in the “Current Personalized Constants” textboxes.



**Automatic Constants Update** - The “Automatic Constants Update” checkbox will enable or disable the Automatic Constants Update feature. This feature automatically calculates and updates the personalized lens constant for the selected doctor and lens based on the records in the table that have “Use Record” selected. The Personalized Lens Constants are also automatically recalculated when new records for the specified doctor and IOL manufacturer and model are created and the “Add to Personalized Lens Table” checkbox on the “Patient Data” Screen is checked.

The number of records used to calculate the constants are also counted and stored automatically upon each recalculation. Automatically Calculated Personalized Constants are only available for use once the number of selected records exceeds the “Min Cases Before Using Personalized Constants” set in Doctor Preferences.

**Min Cases Before Using Personalized Constants** – This value is displayed on the Personalized Constants table. The minimum number of cases used in calculating personalized constants is set in Doctor Preferences. It cannot be edited on this screen.

**Selected Cases** – This value is a counter that notifies the user how many records have been selected in the Use Record column. It cannot be edited on this screen, but will update automatically as records are selected or deselected for use in personalized constant calculations.

## 12.6.2 Procedure

Personalized Lens Constants can be manually calculated by the following procedure:

1. Select the **Patient Data** Screen and either enter a new patient or select an existing patient. Fill out all of the fields in the Patient Info section of the page as described in Section 13.

---

**Note:** Ensure that the nc set in System Preferences matches the corneal refractive index of the Keratometer used to measure the K values.

---

2. Select the **A-Scan** Screen and perform either an Immersion or Contact A-Scan Measurement as described in Sections 14.9 and 14.10.

---

**Note:** DGH does not recommend mixing Contact and Immersion measurement results when calculating personalized lens constants.

---

3. Select the **IOL Calculator** Screen and perform IOL Calculations as described in Section 16. Save the Pre-Operative Data used to perform the calculations by selecting the **Save** button.

4. Perform the lens replacement surgery. Record the actual lens model and power that was used.
5. Enter the Stabilized **Post-Operative** results in the Patient Data Screen. Review the Pre-Operative data and Post-Operative Results and decide if the Procedure should be included in the Personalized Lens Constant Calculation. If the Procedure is a good candidate, select the **Add to Personalized Constants** checkbox.

---

**Note:** Selection of consistent results is the most critical step in Calculating Personalized Lens Constants. Personalized Lens Constants only compensate for consistent errors in predicting the power of the implanted IOL. For this reason, all steps used to measure the patient, calculate the IOL power and perform the lens replacement surgery should be performed as identically as possible. Ultimately it is up to the doctor to decide which results to use, but in general the following guidelines should be observed:

- Do not include results where complications occurred during the procedure.
- Only include results where the same surgical technique was used.
- Do not include results from patients that have undergone refractive eye surgery.
- Only include results where the A-Scan measurement was performed using the same technique, measurement mode (immersion or contact) and instrument.
- Only include results where the Corneal Power measurements were made with the same instrument.

- 
6. After a large population of results has been acquired for a given doctor and IOL model, navigate to the Personalized Lens Constants Table to view the Procedures that were performed by the selected doctor using the selected lens. Review the Procedures that have **Use Record** selected and modify as necessary. See step 5 for considerations on filtering which results to use.
  7. The **New Personalized Constants** group box shows the results of the Personalized Lens Constants calculation based on the Procedures that are currently selected in the table. The “Current Personalized Constants” group box displays the lens constants that are currently being used for the selected Doctor and IOL.
  8. Use the **Up Arrow** to replace the **Current Personalized Constants** with the **New Personalized Constants**. These new Personalized Lens Constants will now be available for use in future IOL power calculations.

Automatically calculating lens constants is similar to the procedure described above except that the lens constants are recalculated as soon as **Add to Personalized Constants Table** checkbox is selected on the Patient Data Screen. See step 5 for more details.

 **WARNING**

Results for Patients That Have Undergone Corneal Refractive Surgery Should Not Be Used To Calculate Personalized Lens Constants.

 **WARNING**

Only A-Scan Measurements Performed Using the DGH 6000 Should Be Used to Calculate Personalized Lens Constants.

 **WARNING**

Only Stabilized Post-Operative Results Should Be Used to Calculate Personalized Lens Constants.

 **WARNING**

Corneal Power Measurements Used For Calculating Personalized Lens Constants Should Be Performed Using the Same Keratometer.

## 13 The Patient Data Screen


Upon startup, the Scanmate application will automatically open to the Patient Data Screen in Search Mode. The Patient Data Screen allows the user to Search, Create, Review and Edit patient records.

Last Name	First Name	ID#	DOB	Gender
Archer	Johnathan	244-9012	10/24/1966	Male
Dahl	Jane	002-0123		Female
Default	Default	Default		Undefined
Doe	John	100-0100	1/9/1954	Male
Smith	Phillip	02/23/2017 16:58:28	12/13/1968	Male
Test Block	Test Block	Test Block		Undefined
Wilkins	Michelle	02/23/2017 16:58:12		Female


### 13.1 Patient Data Screen Controls

The Patient Data Screen operates in three distinct modes of operation: “Search and View”, “Edit / Save Data”, and “New Patient.” The action buttons on the top right-hand side of the screen change availability depending on the current mode of the Patient Data screen. The patient data screen buttons include:


- Open File
- New
- Search
- Edit / Save
- Undo




The **New** button is used to enter new patients in the DGH-Scanmate database. When clicked, the “Last Name”, “First Name”, and “ID#” combo boxes change to plain text boxes and allow the user to enter a new patient.



The **Search** button is used to search for patients in the DGH-Scanmate database. When selected, the “Patient Info” controls are closed for editing and all combo boxes become auto-suggest fields.




The **Edit** button is used to edit patient records in the DGH-Scanmate database. When selected, the “Patient Info” and “Pre-Operative” controls become enabled for editing. The “Post-Operative” controls become enabled for editing only if valid “Pre-Operative” data exists for the selected procedure.



The **Save** button saves New or Modified patient records to the DGH-Scanmate database.



The **Undo** button allows the user to undo the last changes that were made to the Patient’s Record. All fields will revert to the most recent saved state.



The **Open File** button allows the user to open a window to review all records associated with the selected patient.

## 13.2 Entering a New Patient

Select the “New” button and enter the following required fields:

- Last Name
- First Name
- ID Number

The ID Number must be unique for each patient. The Scanmate Software will prevent a new patient from being saved to the database if a patient with the same ID Number already exists. A default unique ID number is created based on the date and time the new patient is entered; this can be replaced with any other numbering system desired. The Patient’s Last Name, First Name and ID Number are required fields and must be entered before the record can be saved.

The following optional fields can also be entered at this time:

- Patient Date of Birth
- Patient Gender
- Doctor
- Comments
- K1 and K2 (K1 and K2 values for OS and OD may be entered from the IOL Calculator Screen and are only needed for IOL Calculations)

- Desired Refraction (Desired Refraction for OS and OD may be entered from the IOL Calculator Screen and are only needed for IOL Calculations)

---

**Note:** It is important that the corneal refractive index (nc) used by the keratometer to calculate the corneal power matches the default corneal refractive index (nc) set in **System Preferences**. The default nc should be set before K1 and K2 are entered for a patient. The nc associated with each set of corneal power measurements can be reviewed or adjusted on the **IOL Calculator Screen**. Failure to adjust the default nc to match the nc of the keratometer used to measure corneal power can result in errors in the IOL Calculation. See Section 16.20 of this guide for more details.

---

Select the **Save** button once all of the desired fields have been entered. A new Procedure will be created once the patient record has been saved. Procedures are used by the Scanmate software to associate a set of patient measurements and calculations (Pre-Operative Data) to a set of surgical outcomes (Post-Operative Data).

### 13.3 Searching for a Patient

Select the **Search** button to search for a patient record that has been saved in the database. Patients can be searched for by Last Name, First Name, or by ID Number.

To search by Last Name, First Name or ID Number, begin typing in the field you wish to search. The software will automatically update the patient table to show all the results that match what has been typed. The desired patient can then be selected from this list.

Once a patient record has been loaded, the Procedures associated with that record can be searched by selecting the Procedure list box.

### 13.4 Editing a Patient Record

Select the desired patient record to be edited following the steps described in Section 13.3. Once the patient record has been loaded, select the **Edit** button. All of the editable fields will change from read-only fields to white, editable text boxes. Once the desired changes have been made, select the **Save** button to save the changed record.

---


**Note:** Once a patient record has been modified, all associated Measurement Files, Video Files and IOL Calculations will be updated to reflect the changes.

---

### 13.5 Creating a New Procedure

Procedures are used by the Scanmate software to associate a set of patient A-Scan measurements and calculations (Pre-Operative Data) to a set of surgical outcomes (Post-Operative Data). Patient records can contain multiple procedures.

To create a new Procedure, select the desired patient record following the steps described in Section 13.3. Once the patient record has been loaded, select “New Procedure” from the Procedure List Box. The pre-operative and post-operative data fields (located on the IOL Calculator Screen) will be blank and ready for the new procedure information to be entered. (Note: To enter data on the Patient Database screen, you must first select the “Edit” button.)



To return to an existing procedure, click on the Procedure List Box and select the desired procedure. Once selected, all of the pre-operative and post-operative data associated with that procedure will be visible.

### 13.6 Entering Pre-Operative Data

Pre-Operative Data is the biometric measurement and procedure data that was used as the basis to perform IOL power calculations for a patient. After a patient has been selected, and a procedure has been created, the Pre-Operative Data is entered on the IOL Calculator Page. The following data must be entered for the procedure:

- K1
- K2
- Desired Rx

The corneal refractive index (nc) associated with the K1 and K2 values can be viewed and changed on the IOL Calculator screen; if it is not changed, the default nc from the System Preferences menu will be used.

### 13.7 Entering Post-Operative Data

The Post-Op Section of the IOL Calculator Screen allows the user to input the stabilized surgical results for the current patient and procedure. This section of the IOL Calculator Screen is only editable if valid Pre-Operative Data has been entered for the selected eye. The following data fields are available:

- **Surgery Date** The date that the lens replacement surgery was performed.
- **Doctor** The surgeon that performed the lens replacement surgery. Personalized constants are performed for this doctor.
- **IOL Mfg** The manufacturer of the implanted IOL

- **Model** The model of the implanted IOL
- **Power** The power of the implanted IOL
- **Sph** The Stabilized Post-Operative Spherical Refractive Error
- **Cyl** The Stabilized Post-Operative Cylindrical Refractive Error
- **Add to Personalized Constants Table**  Add to Personalized Constants Table?  
 Selecting this checkbox indicates that the Doctor would like to include these Results in the Calculation of Personalized Lens Constants. See Section 12.6 for more information about Personalized Lens Constants.



**WARNING**

Results for Patients That Have Undergone Corneal Refractive Surgery Should Not Be Used To Calculate Personalized Lens Constants.



**WARNING**

Only A-Scan Measurements Performed Using the DGH 6000 Should Be Used to Calculate Personalized Lens Constants.



**WARNING**

Only Stabilized Post-Operative Results Should Be Used to Calculate Personalized Lens Constants.




**WARNING**

Corneal Power Measurements Used For Calculating Personalized Lens Constants Should Be Performed Using the Same Keratometer.

## 14 The A-Scan Screen

The A-Scan Screen allows the user to perform and review A-Scan Measurements for the currently selected patient. Refer to Section 13 for details on how to load or create patient records.



The default settings for this page are dependent on the currently selected Doctor and Operator. The default settings can be changed by selecting **Preferences** → **Operator** or **Preferences** → **Doctor**.

### 14.1 Selecting the Operator

The name of the Operator performing the measurement can be selected using the Operator drop down box. If the operator is not listed in the drop-down box, add the new Operator by selecting **Preferences** → **Operator** and clicking the **New** button.

### 14.2 Selecting the Lens and Vitreous Type

The lens and vitreous type of the patient can be selected using the drop-down box and radio buttons directly under the operator selection box. The currently selected lens and vitreous type is displayed in the top left-hand corner of the waveform display screen.

### 14.3 Selecting OD or OS

The eye currently being measured (OD or OS) can be selected by pressing the **Change Eye** button. The currently selected eye is shown in the top left-hand corner of the waveform display screen.



### 14.4 Measurement Modes

The DGH Scanmate software can be configured for Automatic or Manual measurement mode. The Automatic and Manual Measurement modes are selected using the radio buttons to the right of the waveform display. Scans performed in Automatic Measurement Mode can be replayed in Manual Measurement Mode, allowing the user to adjust measurement points. Likewise, scans performed in Manual Measurement Mode can be replayed in Automatic Measurement Mode, allowing the automatic entry of measurements into the measurement bank.



### 14.5 Automatic Measurement Mode


In Automatic Measurement Mode, the DGH 6000 uses a pattern recognition algorithm to determine when the probe is properly aligned. The alignment algorithm first looks for the proper pattern of reflectance peaks from each expected interface. Then the algorithm examines the retinal spike for particular signal characteristics which are produced only by retinal interfaces. Waveforms that meet the alignment criteria are ranked and put into the measurement bank. Alignment is ranked using a “star” system. As measurements enter the measurement bank, they are given a 1-star, 2-star or 3-star rank, with 3-star being the best. 3-star ranked measurements will replace those in the bank with a lower ranking.

Fr#	AXL	ACD	LT	VCD	Rank
6	25.68	3.74	3.77	18.17	***
7	25.68	3.74	3.77	18.17	***
8	25.65	3.71	3.76	18.18	***
9	25.62	3.69	3.75	18.18	***
10	25.58	3.65	3.75	18.18	***
11	25.55	3.63	3.73	18.19	***
12	25.58	3.66	3.72	18.20	***
13	25.60	3.68	3.72	18.20	***


In order to meet the measurement criteria, eight measurements must be obtained that meet the ranking (3-Star) and standard deviation requirements. Standard deviation requirements can be viewed and adjusted by selecting **Preferences → Doctor**.

Once the measurement bank is filled with 3-star measurements, and all measurements are within the standard deviation parameters, measurement will automatically stop. The entries in the bank are used to calculate the measurement average. Refer to Sections 14.9 and 14.10 for instructions on performing an A-Scan.

If 3-star alignment is not achieved within the alignment timeout period, a warning message will appear stating that alignment accuracy requirements will be reduced. 1-star, 2-star and 3-star measurements will be accepted for calculating the average.



If alignment continues to be undetected, a warning message will appear stating that alignment accuracy requirements will be disabled, placing the system in “Estimate Mode”. All measurements, regardless of alignment, will be accepted for calculating the average.



If measurements are obtained that do not show signs of alignment, they will be given an “EST” ranking.

168	23.72	3.57	3.96	16.19	EST
-----	-------	------	------	-------	-----

---

**Note:** Measurements with an “EST” ranking do not show signs of alignment. Star ranked measurements should be given preference over “EST” ranked measurements when possible.


---

If eight valid measurements are not obtained prior to the video buffer expiring, the scan will time out. Two long beeps will signal that the scan has timed out.

#### 14.6 Manual Measurement Mode

In Manual Measurement Mode, scans must be stopped by the operator and measurements must be manually selected and entered into the measurement bank. Manual Mode scans are only stopped when the user selects to stop the scan or if the video buffer expires. Refer to Sections 14.9 and 14.10 for instructions on performing an A-Scan.

After an A-Scan video has been acquired, it can be reviewed in Manual Measurement Mode by using the ← and → cursor keys or moving the slider below the waveform display to scroll through the captured waveforms.



On the desired frame(s), the measurement can be added to the bank by pressing the **Enter** key (on the keyboard), **Right Foot Pedal** (optional accessory) or by right-clicking on the measurement bank or waveform display and selecting **Insert Measurement**.

Measurements that are manually inserted are given a “MAN” rank. Measurements can be removed from the measurement bank by right-clicking on the measurement with the mouse and selecting Remove Measurement.

110	23.78	3.70	3.71	16.37	Man
-----	-------	------	------	-------	-----

---

**Note:** Measurements with an “MAN” ranking are not evaluated for alignment. It is the user’s responsibility to ensure that measurements acquired in this manner are properly aligned.

---



On a desired frame, the measurement points can be set by the user. In Contact Mode three red arrows represent measurement points that can be selected; they represent the anterior and posterior lens surfaces (L1 and L2) and the retina (R). In Immersion Mode, the corneal peak is also visible in the waveform and is indicated by an additional arrow at the leftmost echo.

The red arrows can be moved horizontally by clicking on an arrow with the mouse, holding the mouse button, and moving the mouse until the arrow is in the desired position. When a measurement arrow is released the arrow snaps to the closest intersection of the measurement level and the displayed waveform. The measurement points can only be adjusted when the video is paused.

The horizontal green line represents the measurement level. In manual mode, it can be adjusted up and down while the video is paused or playing.

The three sliders control the gain at the three measurement points (L1, L2, and R). They can be adjusted independently while the video is paused or playing. The amount of gain is displayed as a percentage of the incoming waveform amplitude.

## 14.7 Selecting Immersion or Contact Measurement Modes

The DGH 6000 is capable of performing both “direct contact” and “water immersion” measurements. Direct contact means that the ultrasonic probe must be applanated to the cornea. This method is the quickest and most convenient way of taking measurements.

Because there is the possibility of indenting the cornea using the Contact method, the measured axial length may be shorter than its actual value. To improve accuracy, the DGH 6000 has a pattern recognition program that rejects any measurement which contains a significant amount of corneal indentation; however, some degree of indentation may still exist. The skill and technique of the operator will reduce possible corneal indentation. The Compression Lockout Sensitivity setting determines how much indentation will be tolerated (see Section 14.8).

When using the “Immersion” mode, measurements are obtained in a column of water while the probe is suspended above the eye. Since the probe does not actually touch the eye, corneal compression is not an issue. Immersion mode increases both measurement accuracy and precision. To completely eliminate corneal indentation, “Immersion” mode should be used.


The Contact and Immersion Measurement modes are selected using the radio buttons to the right of the waveform display screen.



When performing measurements in immersion mode, the waveform will appear as a yellow line.



When performing measurements in contact mode, the waveform will appear as a white line.



## 14.8 Setting the Corneal Compression Lockout Sensitivity

When in Contact Measurement Mode, the Operator can control the Corneal Compression Lockout Sensitivity.

The Comp Sens slider to the right of the waveform display screen adjusts the amount



of corneal compression that the software will allow when taking measurements in contact mode. Setting this slider higher will allow less compression, making measurements more accurate but also more difficult to obtain. The optimal setting for Comp Sens will be different for each probe and will depend on the operator’s technique.

The Corneal Compression Lockout Sensitivity can be adjusted before replaying videos or on any paused frame to see the impact on the measurements. A red/green indicator next to the Comp Sens slider will show whether a measurement can be taken from a frame: if the indicator is red, there is too much compression for a measurement at that sensitivity level. If the indicator is green, a measurement can be taken. When properly set, the Comp Sens control should limit corneal compression in Contact Mode to 0.15 mm or less.


## 14.9 Performing a Contact Measurement

Before each biometry procedure, the ultrasound probe must be cleaned. (See Section 20.2 for Cleaning and Disinfection Instructions.)

---

**Note:** † Indicates a command that may be executed through A-Scan voice controls. Refer to Section 14.14 for more information.


---

1. Launch the Scanmate Software Application as described in Section 11.
2. Select the **Patient Data** Page and either enter a new patient or select an existing patient. Fill out all of the fields in the Patient Info section of the page as described in Section 13.
3. Navigate to the **A-Scan** page and select the eye to be measured (OD or OS) by clicking the **Change Eye**† button.
 
4. Configure the unit for contact mode by selecting the **Contact** radio button to the right of the waveform display.
 
5. Select the name of the Operator performing the measurement. (If the operator is not listed in the drop-down box, add the new operator by selecting **Preferences** → **Operator** and clicking the **New** button.)
6. Select the correct Lens Type and Vitreous type for the eye to be measured.
7. Adjust the Compression Lockout Sensitivity (Comp) track bar to the desired setting. Increasing the Compression Sensitivity will decrease the amount of error due to corneal compression, but will make it more difficult to obtain measurements. See Section 14.8 for more details.
8. Position the display for easy visibility during patient examination.
9. Seat or recline the patient in a comfortable position for both the Patient and Operator. Use a firm, comfortable head rest to prevent head movement during the exam.
10. Anesthetize the cornea and ask the patient to fixate with their fellow eye on a spot on the ceiling or wall.
11. Press the **Start Scan**† button. Alternatively, clicking the **Space Bar** (on the keyboard) or **Left Foot Pedal** (optional accessory) will also begin a scan. A
 

waveform will be displayed on the screen, but the unit will not begin recording.

12. The hand holding the transducer can be stabilized on the cheek or forehead of the patient. This will help minimize indentation of the cornea or excessive movement.

13. The transducer can now be applanated to the cornea. Applanation should be observed from a vantage point that allows the operator to look across the topography of the eye. The transducer should be perpendicular to the surface of the eye as it approaches the cornea. This method aids in seeing the exact moment when the transducer contacts the cornea so



indentation of the cornea will be minimized. Applanation should be made as close to the visual axis as possible. An aid in this can be to point the back side of the transducer at the fixation target while touching the front portion of the transducer to the center of the cornea.

---

**Note:** When applanating the transducer, be particularly careful that contact is made with the center of the transducer tip.

---

14. At the moment of applanation, the unit will start beeping at a rate of approximately one beep per second. This means that proper coupling has been achieved and sound pulses are being transmitted through the eye. If the beeping stops, then the transducer is no longer applanated to the eye and a slight pressure should be applied to restore applanation. The DGH 6000 will only record the live waveforms and perform measurements when it detects that proper coupling has been achieved.

**\*\*\* IMPORTANT \*\*\***

In order to obtain measurements, the transducer must be in proper alignment. Audible feedback from the unit is a key element for achieving proper alignment. Measurements are best obtained by concentrating on the audible feedback while visually aligning the transducer. It is not recommended that the operator view the waveforms on the graphic display during the actual measurement cycle. The interpretation of the audible tones that are emitted during the measurement cycle is explained below.

15. The alignment of the transducer may need slight adjustments; keep the tip in contact as the angle is changed slightly. As the transducer is brought closer to alignment, the beep rate will increase.

16. When the transducer is properly aligned, a beep with a higher tone will sound. The axial length measurement can be observed on the graphic display along with its echo waveform. Measurements for the anterior chamber depth and lens thickness are also given. If the unit is configured for the “Aphakic” lens type, only the axial length will be given.

In **Auto Measurement Mode**, each high-pitched beep indicates a successful measurement being recorded in the measurement bank. Measurements automatically placed in the bank are given a 1, 2 or 3-star alignment ranking. Scans can also be aborted manually as described below.

In **Manual Measurement Mode**, a measurement is only recorded in the measurement bank if the user presses the **Enter** key (on the keyboard) or **Right Foot Pedal** (optional accessory). Measurements manually placed in the bank are given a “MAN” ranking.

17. Scanning is stopped as follows:

In **Auto Measurement Mode**, the scan will stop automatically as soon as eight (8) measurements are obtained that meet the alignment ranking and standard deviation requirements.

In **Manual Measurement Mode**, the user must select the **Stop Scan†** button to stop the scan. Alternatively, clicking the **Space Bar** (on the keyboard) or **Left Foot Pedal** (optional accessory) will also stop a scan.

18. If a measurement is not obtained successfully when applanation occurs and the beep rate is still approximately once per second, remove the transducer from the cornea and position the transducer at a slightly different location or angle. Repeat this removal and repositioning procedure until either a successful measurement or a very rapid beep rate is obtained.

19. At the higher beep rate, measurements can be obtained by slightly tilting the back end of the transducer while the tip remains stationary and applanated at the same location. As soon as proper alignment is achieved, a measurement will be obtained provided that there is no corneal indentation present.

20. If a lower tone beep is heard during the alignment procedure, cornea indentation may be occurring. The DGH 6000 uses a special algorithm that prohibits any measurements from being obtained if the cornea is significantly indented. The Compression Lockout Sensitivity (Comp) track bar adjusts the sensitivity of the unit in detecting corneal compression. If low-pitched beeps are heard, maintain the transducer alignment and slowly withdraw the transducer from the cornea until measurements are obtained.



21. If 8 valid measurements have not been obtained before the video buffer is full, the scan will time out as indicated by two long beeps. A message will appear

indicating that the measurements did not meet the measurement selection criteria or the standard deviation requirements.

---

**Note:** The default video buffer duration is 1024 frames. This can be adjusted in the **Preferences** → **System** menu.

---

22. All of the waveforms captured during the exam can be reviewed on the display by pressing the → or ← cursor keys. If the currently selected waveform is associated with a measurement in the measurement bank, the corresponding record in the measurement bank will be highlighted.

23. The video buffer can be played or rewound to the beginning by clicking on the arrow buttons in the bottom left corner. The slider bar beneath the waveform display can also be used to select an image for viewing from the video buffer series.

24. If a measurement is observed in the measurement bank which is suspect because it does not agree with the other measurements, it may be removed or replaced with a different measurement:

a) Select the Manual Mode radio button (if not already in Manual Measurement Mode).




b) Using the mouse, **RIGHT-CLICK** the measurement in the measurement bank you wish to remove and select “Remove Measurement” from the menu.

c) Use the → and ← cursor keys to scroll through the captured waveforms and add the desired measurement to the measurement bank by **RIGHT-CLICKING** on the measurement bank and selecting “Insert Measurement” from the menu.

d) Select the Auto Measure radio button if it is desired to continue in Automatic Measurement Mode.

25. After obtaining the desired measurements for the first eye (OD or OS), the operator may obtain measurements on the same patient’s other eye, if desired, by clicking on the **Change Eye**† button and repeating the steps above.




---

**Note:** The measurement data for the first eye will remain stored in memory while measurements are being obtained on the second eye.

---

26. After measurements have been taken, the operator may save a video that stores all of the waveforms captured during the exam by selecting **File** → **Save** → **A-Scan Video**† from the top menu bar. Alternatively, a measurement file that only stores the eight (8) entries in the measurement bank and their associated waveforms can be saved by clicking the **Save**† button or by selecting **File** → **Save** → **A-Scan Measurements**.



27. The video buffer and associated measurements will remain in memory for each eye (OS or OD) until either a new scan is performed  for that eye, or a new patient is selected. The **Rewind** button clears the measurement bank of the currently selected eye, but does not clear the video buffer.

## 14.10 Performing an Immersion Measurement

---

**Note:** The Prager Shell included with the DGH 6000 comes with instructions for use and images illustrating the catheter tubing connection. Please refer to these sheets for more detailed instructions regarding the immersion shell.

---



Before each biometry procedure, both the ultrasound transducer and the immersion shell must be cleaned. This is in addition to disinfection procedures that should be performed between patients. (See Sections 20.2 and 20.3 for Cleaning and Disinfection Instructions.) The asepsis procedure most commonly employed with a biometry transducer is a wipe-down with an alcohol prep. The immersion shell can be cleaned by removing the catheter tubing and washing using a bactericidal liquid soap followed by a thorough hot water flush; some practitioners insure asepsis by including a soak in hydrogen peroxide, alcohol, or Cidex. It is recommended that busy biometry departments should acquire two or three shells, allowing one or more shells to be in the cleaning cycle while another one is in use.

DGH recommends using fresh saline for each patient. The catheter tubing provided with the immersion shell has a backflow check valve, so it can be attached to either a syringe or a sterile BSS squeeze bottle. If using tubing that does not have a check valve, only a syringe may be used.

---

**Note:** † Indicates a command that may be executed through A-Scan voice controls. Refer to Section 14.14 for more information.

---

1. Launch the Scanmate Software Application as described in Section 11.
2. Select the **Patient Data** Page and either enter a new patient or select an existing patient. Fill out all of the fields in the Patient Info section of the page as described in Section 13.
3. Navigate to the **A-Scan** page and select the eye to be measured (OD or OS) by pressing the **Change Eye**† button. 
4. Configure the unit for immersion mode by selecting **Immersion** radio button to the right of the waveform display. 
5. Select the name of the Operator performing the measurement. (If the operator is not listed in the drop-down box, add the new operator by selecting **Preferences** → **Operator** and clicking the **New** button.)

6. Select the correct Lens Type and Vitreous type for the eye to be measured.
7. Position the unit for easy visibility during patient examination.
8. Seat or recline the patient in a comfortable position so that the head can be positioned in a nearly horizontal plane. Use a firm, comfortable head rest to prevent unwanted head movement.
9. Insert the biometry probe into the clean and disinfected shell. Gently insert the transducer into the shell until it reaches the auto-stop point. Verify that the transducer tip is even with the line scored on the barrel of the shell. Once the correct position has been found, tighten the white nylon set screw to hold the transducer in place. Do not overtighten the set screw. The software will work only if the probe is in the proper position.
10. Connect the catheter tubing to the immersion shell filler port and to the saline reservoir (either a syringe or bottle).
11. Administer routine topical anesthesia to the patient's eye.
12. Place a towel on the patient's shoulder and rest the saline reservoir on the towel. Hold onto the transducer and shell in preparation for the insertion, with the filler port and transducer cord oriented in a position that is comfortable for the operator and the patient.
13. Press **Start Scan**† button. Alternatively, clicking the **Space Bar** (on the keyboard) or **Left Foot Pedal** (optional accessory) will also begin a scan.
 


**Start Scan** • 1 ) )
14. Direct the patient to look upward. Pull the patient's lower eye lid down and insert the flared rim inside the lid (the lower portion of the shell will make contact with the sclera while the upper part of the shell will be held away from the eye). Then direct the patient to look straight ahead with the other, uncovered eye, toward a fixation point. Lift the patient's upper eyelid and gently pivot the upper portion of the shell into the upper fornix, making sure by close inspection that it is in the fornix and not atop a fold in the conjunctiva. This pivotal motion avoids contact with the cornea and insures centration of the device around the limbus.
15. Rest the left palm holding the immersion shell on the patient's forehead to reduce shell pressure on the eye. In this position, the palm acts as the fulcrum or pivot point for the shell. If necessary, the operator can stabilize the shell with the other hand and make micro-adjustments. A tissue may be placed on the temporal canthus to catch any excess saline.

16. Pick up the saline reservoir from its place on the patient's shoulder and slowly inject the saline into the shell. As soon as the liquid fills the shell sufficiently to reach the tip of the transducer (about 2cc), a beep rate of approximately once per second will be heard. This means that proper coupling has been achieved and sound pulses are being transmitted through the eye.

**\*\*\* IMPORTANT \*\*\***

In order to obtain measurements, the transducer must be in proper alignment. Audible feedback from the unit is a key element for achieving proper alignment. Measurements are best obtained by concentrating on the audible feedback while viewing the transducer with shell as it is aligned with the visual axis. It is not recommended that the operator view the waveforms on the graphic display during the actual measurement cycle. The significance of the audible tones that are emitted during the measurement cycle is explained below.

17. If the transducer is in proper alignment when the saline is injected, measurements will be obtained almost instantaneously as indicated by a beep with a higher tone. The axial length measurement can be observed on the graphic display along with its echo waveform. Measurements for the anterior chamber depth and lens thickness are also given. If the unit is configured for the "Aphakic" lens type, only the axial length will be given.



18. As soon as proper alignment is achieved, a measurement will be obtained as indicated by a beep with a higher tone. Each high-pitched beep indicates a measurement being recorded in the measurement bank.

**In Auto Measurement Mode**, each high-pitched beep indicates a successful measurement being recorded in the measurement bank. Measurements automatically placed in the bank are given a 1, 2 or 3-star alignment ranking.

**In Manual Measurement Mode**, a measurement is only recorded in the measurement bank if the user presses the **Enter** key (on the keyboard) or **Right Foot Pedal** (optional accessory). Measurements manually placed in the bank are given a "MAN" ranking.

19. Scanning is stopped as follows:

In **Auto Measurement Mode**, the scan will stop automatically as soon as eight (8) measurements are obtained that meet the alignment ranking and standard deviation requirements. Scans can also be aborted manually as described below.

In **Manual Measurement Mode**, the user must select the **Stop Scan†** button to stop the scan. Alternatively, clicking the **Space Bar** (on the keyboard) or **Left Foot Pedal** (optional accessory) will also stop a scan.

20. If a measurement is not obtained instantaneously and the beep rate remains approximately once per second, then the transducer is out of alignment. Move the shell around the sphere of the eye to bring the line of the transducer closer to the visual axis. The beep rate will increase as the transducer approaches the visual axis.
21. If eight (8) valid measurements are not obtained prior to the video buffer expiring, the scan will time out as indicated by two long beeps. A message will appear indicating that the measurements did not meet the measurement selection criteria or the standard deviation requirements.


---

**Note:** The default video buffer duration is 1024 frames. This can be adjusted in the **Preferences → System** menu.

---


22. If the unit has trouble obtaining measurements, then the transducer alignment procedure described in steps 17-20 should be repeated until successful. If transducer alignment is not successful, the shell should be removed from the patient's eye (see step 23) and the procedure should be repeated from step 13.
23. After measurements are obtained, remove the shell from the patient's eye. Direct the patient to look straight ahead during removal. Raise the patient's upper eyelid, which releases the top part of the shell from under the eyelid. Then pivot the shell downward, directing the patient to continue to look straight ahead. Then pull the shell away from the eye without contacting the cornea. Upon the initial release, the remaining contents of the shell (1-3cc of liquid) will spill down the patient's cheek. Be prepared with a towel or tissue.
24. All of the waveforms captured during the exam can be reviewed on the display by pressing the → or ← cursor keys. If the currently selected waveform is associated with a measurement in the measurement bank, the corresponding record in the measurement bank will be highlighted.
25. The video buffer can be played or rewound to the beginning by clicking on the arrow buttons in the bottom left corner. The slider bar beneath the waveform display can also be used to select an image for viewing from the video buffer series.
26. If a measurement is observed which is suspect because it does not agree with the other measurements, the measurement may be removed and replaced as follows:

- a) Select the Manual Mode radio button (if not already in Manual Measurement Mode).



- b) Using the mouse, **RIGHT-CLICK** the measurement in the measurement bank you wish to remove and select “Remove Measurement” from the menu.
- c) Use the → and ← cursor keys to scroll through the captured waveforms and add the desired measurement to the measurement bank by **RIGHT-CLICKING** on the measurement bank and selecting “Insert Measurement” from the menu.
- d) Select the Auto Measure radio button if it is desired to continue in Automatic Measurement Mode.

27. After obtaining the desired measurements for the first eye (OD or OS), the operator may scan the same patient’s other eye by clicking the “Change Eye” button and repeating the steps above.




---

**Note:** The measurement data for the first eye will remain stored in memory while measurements are being obtained on the second eye.

---

28. After measurements have been taken, the operator may save a video that stores all of the waveforms captured during the exam by selecting **File → Save → A-Scan Video†** from the top menu bar. Alternatively, a measurement file that only stores the eight (8) entries in the measurement bank and their associated waveforms can be saved by clicking the **Save†** button or by selecting **File → Save → A-Scan Measurements**.




29. The video buffer and associated measurements will remain in memory for each eye (OS or OD) until either a new scan is performed for that eye, or a new patient is selected. The **Rewind** button clears the measurement bank of the currently selected eye, but does not clear the video buffer.



### 14.11 Saving a Measurement File


A Measurement File is a group of up to eight (8) A-Scan measurements and their associated waveforms for both OS and OD. Measurement Files can be saved by selecting the “Save” Button on the upper right-hand corner of the A-Scan Screen. The **Save** Button will also save all changes that have been made to the Patient Data or IOL Calculator since the last save.



Alternatively, Measurement Files can be saved using the toolbar at the top of the screen by selecting **File → Save → A-Scan Measurements**.


### 14.12 Saving a Video Buffer

From the toolbar at the top of the screen, select **File → Save → A-Scan Video** to save a video recording of the current exam.



### 14.13 Reviewing Measurements and Videos

From the toolbar at the top of the screen, select **File → Open** or click the **Open File** button to review previous exam information captured for the current patient. A window will open showing saved measurement files and videos for the currently loaded patient. Double-clicking on the desired measurement file or video will open the record for review. Videos can be played back in the measurement window and the AXL, ACD and LT measurements can be recalculated.



### 14.14 Voice Controls

Some controls can be operated via verbal commands from the user. Voice control can be enabled/disabled through the **Preferences → System** menu. Voice control can be configured to begin upon pressing the “Enter” key or by a verbal command of “Hey Flex”. Voice control sensitivity can be adjusted in this menu as well. The following verbal commands are recognized by the Scanmate software while performing an A-Scan.


User Verbal Commands	Software Audible Response	Software Command Executed
“Hey Flex”	Tone: Double-Beep	Voice control activated. The software is ready to receive verbal commands.
“Start Scan”	Verbal: “Starting Scan”	The A-scan will start.
“Stop Scan” *	Verbal: “Scan Stopped”	The A-scan will stop.
“Change Eye”	Verbal: “Eye changed to OD/OS.”	The eye being scanned is changed. The audible response will state if OD or OS is being scanned.

“Save Measurements”	Verbal: “Measurements Saved”	The currently displayed frame is saved.
“Save Video”	Verbal: “Video Saved”	The currently displayed video is saved.
“Open File”	N/A	A window is opened showing A-Scan measurements and videos for the currently loaded patient.
“Yes” *	N/A	Selects ‘Yes’ for a dialog box prompts.
“No” *	N/A	Selects ‘No’ for a dialog box prompts.

\* Indicates a command that DOES NOT require “Hey Flex” or an “Enter” key press to be initiated.

## 15 The B-Scan Screen

The B-Scan Screen allows the user to perform and review B-Scan exams for the currently selected patient. Please refer to the DGH 8000 Scanmate-B or Scanmate Flex User Guide for more information on using the B-Scan Screen.



## 16 The IOL Calculator Screen

The IOL Calculator Screen allows the user to perform and review intraocular lens power calculations for the currently selected patient. If available, measurements and the selected lens and vitreous type are automatically loaded from the A-Scan page. Information from the Patient Data screen, including Patient Name, Patient ID, and Doctor are also automatically loaded.

The screenshot displays the IOL Calculator interface for a patient named Jane Doe, dated 06/29/2017 at 11:16:43. The interface is split into two main sections: OD (Right Eye) and OS (Left Eye). Each section includes a 'Pre-OP' and 'Post-OP' toggle, scan date, scan type, source, and total number of scans. Below these are tables for measurements (Meas), standard deviations (SD), and refractive indices (nc). The OD eye has measurements: ACD 3.74 mm, LT 3.67 mm, AXL 25.54 mm, K1 (Flat) 40.00 D, K2 (Steep) 40.00 D, and Desired Rx 0.00 D. The OS eye has measurements: ACD 3.69 mm, LT 3.75 mm, AXL 25.62 mm, K1 (Flat) 40.00 D, K2 (Steep) 40.00 D, and Desired Rx 0.00 D. At the bottom of each section is a table for IOL selection, showing IOL Power and Refraction for three IOLs. The OD eye table shows IOL Power values of 16.00, 16.50, 17.00, 17.50, 18.00, 18.50, and 19.00, with corresponding Refractions of +1.32, +0.92, +0.52, +0.12, -0.28, -0.68, and -1.08. The OS eye table shows IOL Power values of 16.00, 16.50, 17.00, 17.50, 18.00, 18.50, and 19.00, with corresponding Refractions of +1.16, +0.76, +0.36, -0.04, -0.44, -0.84, and -1.24. The A-Const is set to 118.00 for both eyes.

The Pre-Op section of the IOL calculator contains fields for all of the necessary patient data inputs for performing an IOL calculation. Any inputs for performing an IOL Power Calculation that are either missing or out of range are automatically highlighted in red. The results of the IOL calculation (Target IOL Power, Emmetropic Power, and Range of IOL Power and Refractions) are displayed for up to three (3) IOLs for each eye.

## 16.1 IOL Calculator Screen Controls


Depending on the selections made by the user, up to four (4) action buttons are available at the top right corner of the IOL Calculator Screen.



The **K-Corr** button launches the Post Refractive Surgery Calculator. The “Post Refractive” checkbox must be selected for at least one eye (OS or OD) for this button to be available.



The **Save** button saves the pre-operative data for the current patient and current procedure. When clicked, the pre-operative data on the IOL Calculator screen is sent to the patient data screen and will become part of that patient file. In addition, measurement data on the A-Scan Page and patient data is automatically saved and the IOL information is automatically saved with the pre-op data.



The **Report** button will create an IOL Calculator Report. You can also use the menu bar to go to **Reports → IOL Calculator Report**. When clicked, a window will open displaying the report. The report can then be saved to the database, exported as a PDF file or sent to a printer.

## 16.2 Patient Information

The Name and ID Number of the patient currently selected on the Patient Data Screen are displayed in the top left corner of the IOL Calculator Screen. This information can be verified to confirm that the calculation is being performed on the correct patient.

## 16.3 Selecting the Doctor


The doctor performing the IOL Calculations can be selected using the **Doctor** drop-down box. This field will default to the Doctor that was selected on the Patient Data Screen. If the doctor is not listed in the drop-down box, add the new doctor by selecting **Preferences → Doctor** and clicking the **New** button. When a doctor is selected, the doctor’s preferred IOLs and Formulas will automatically be loaded into the IOL Calculator. See Section 12.5 for more information on configuring preferred IOLs and formulas.

## 16.4 Eye Selection

The Scanmate software will automatically perform IOL Power Calculations once all necessary parameters for the subject eye(s) have been input.

## 16.5 Post Refractive Selection

Select the **Post Refractive** checkbox if the patient has previously undergone refractive eye surgery. Selecting this checkbox enables the post refractive calculator which assists the doctor in predicting the true corneal power of eyes that have undergone corneal refractive surgery. See Section 16.16 and 16.17 for more details on using the Post Refractive Calculator.



## 16.6 Measurement Source

The **Source** drop-down box allows the user to select the source of the Pre-Operative Measurement values. The following options are available:



- Unit Average
- Unit Single
- User Input

When “Unit Average” is selected, the average and standard deviation of the Anterior Chamber Depth, Lens Thickness, and Axial Length measurements are displayed in the pre-operative data section of the IOL Calculator. The total number of measurements used to create the average is displayed in the “Total #” field.

When “Unit Single” is selected, the Anterior Chamber Depth, Lens Thickness, and Axial Length measurements for the currently selected waveform on the A-Scan page is displayed in the pre-operative data section of the IOL Calculator.

When “User Input” is selected, the user is able to manually enter the desired Anterior Chamber Depth, Lens Thickness, and Axial Length measurements.

The Measurement Source field will default to the current doctor’s default measurement source. The measurement type, lens type and vitreous type will be imported from the A-Scan screen when either “Unit Average” or “Unit Single” is selected. When “User Input” is selected, the lens type and vitreous type fields become drop-down boxes for the user to select and the measurement type field is removed.

## 16.7 Corneal Refractive Index

The Corneal Refractive Index (nc) of the keratometer used to calculate the corneal power (K1 and K2) for the current patient should be entered in this field. This field will automatically be completed with the default corneal refractive index (nc) set in System Preferences, but it can be changed by the user on this screen.

---

**Note:** Failure to adjust the default nc to match the actual nc of the keratometer used to measure corneal power can result in errors in the IOL Calculation. See Section 16.20 of this guide for more details.

---

## 16.8 Corneal Power (K1 and K2)

The corneal power values (K1 and K2) for the current patient are entered in these fields. When the **Save** button is selected on the IOL Calculator Screen, these values will automatically be saved in the “Pre-Op” section of the IOL Calculator Screen for the current patient and procedure. These fields must be completed before performing IOL calculations. Corneal power can be entered in either millimeters (mm) or diopters (D); the software will automatically determine which units to use, so the user only needs to enter the number. K1 and K2 must be entered in the same units.

For Post Refractive calculations, these fields will display the Corrected Corneal Power (Kcorr) calculated by the Post Refractive Surgery Calculator. Kcorr is a **read-only** value, which can only be modified by changing the inputs to the Post Refractive Surgery Calculator. Refer to Sections 16.16 and 16.17 for more information on the Post Refractive Surgery Calculator.

## 16.9 Desired Refraction

The desired surgical refractive result (Desired Rx) is entered in this field. When the Save button is selected on the IOL Calculator Screen, the entered value will automatically be saved in the “Pre-Op” section of the IOL Calculator Screen for the current patient and procedure.

## 16.10 Selecting IOLs

These tabs (IOL 1, IOL 2, and IOL 3) allow the user to select one of three preferred IOLs on which to base the IOL implant power calculation. This selection will default to the preferred IOLs for the currently selected Doctor. If the desired IOL is not one of the three preferred IOLs, it can be selected using the drop-down lists. First select the manufacturer from a list of all the available manufacturers for IOLs entered for the current doctor. The model drop-down list will then offer a choice of all the IOL models entered for that manufacturer. See Section 12.5 for more information on configuring preferred IOLs and formulas.

## 16.11 Selecting IOL Power Calculation Formulas

The formula selection drop-down box allows the user to choose a formula to use to perform an IOL Power Calculation for the selected lens. Each pre-set IOL has a preferred formula associated with it that will appear automatically in this box, but all other formulas can be selected from the drop-down list. If the Post Refractive Calculator has been used to calculate Kcorr, only the appropriate formulas will be available for selection.

## 16.12 Selecting Personalized Lens Constants

The **Personalized** checkbox allows the user to use a personalized lens constant for the selected formula. This checkbox is only selectable if the lens constant for the selected formula can be personalized, and the “minimum number of cases” has been reached using the currently selected IOL. (See Section 12.6 for more information on how to calculate personalized lens constants.)



## 16.13 K Correction Method Selection

This drop-down box selects between different methods of predicting the true corneal power of eyes that have undergone corneal refractive surgery. This field is only available if the Post Refractive Calculator has been used. If the K Correction Method is changed to a method that requires information that has not been completed, the Post Refractive Calculator will reopen to request the additional inputs. If the new method chosen does not require more information, the Kcorr value will automatically be updated to the result of the method.

## 16.14 IOL Calculator Outputs

Once all necessary inputs are entered, the DGH Scanmate will perform IOL calculations using the selected formula for each of the six (6) selected IOLs on the screen. Any necessary inputs for performing an IOL Power Calculation that are either missing or out of range are automatically highlighted in red.


The Target Power is the ideal result. However, since IOLs are only manufactured in discrete power increments, an IOL must be selected having a power that comes closest to the Target Power. To aid the operator in selecting an IOL with the most appropriate power, the DGH Scanmate performs the following calculations:

- An IOL power that predicts Emmetropia.
- A target IOL power that is based on the desired postoperative refraction.
- A selected IOL having a power that comes closest to the target IOL power.
- A postoperative refraction predicted by using the selected IOL power.
- Postoperative refractions predicted by using IOL powers for five power increments above and below the selected IOL power.

The power increments are set for each IOL in the IOL configuration table in Doctor Preferences, by selecting  $\frac{1}{2}$  or  $\frac{1}{4}$  diopter.

## 16.15 Saving Pre-Operative Data

The **Save** button allows the user to save the pre-operative data for the current patient and current procedure. When clicked, the pre-operative data on the IOL Calculator screen is sent to the patient data screen; measurement data on the A-Scan Page is saved if unsaved; and the IOL information is saved with the pre-operative data.



## 16.16 Post Refractive Calculation Methods

It has been well established that the corneal power for patients that have undergone corneal refractive eye surgery cannot be accurately measured using current keratometry or topography methods.

The Post Refractive Calculator included with the DGH Scanmate helps predict the true corneal power for patients that have undergone corneal refractive surgery. The DGH Scanmate is equipped with the following methods for calculating the corrected corneal power (Kcorr):

- History Derived
- Clinically Derived
- Refraction Derived
- Contact Lens Over-Refracton
- Double K

These methods have been implemented as described in H. John Shammas's book "Intraocular Lens Power Calculations" published by Slack Incorporated 2004. The mathematical formulas used for each of these methods can be found in Appendix C of this User Manual.

**Post Refractive K Calculator**

Calculate OD      Method

**OD**

Inputs

Pre Refractive Surgery

nc 1.3375

K1  K2

Spectacle Refraction (S.E.)  D

Post Refractive Surgery

nc 1.3375

K1 7.00 mm K2 7.00 mm

Spectacle Refraction (S.E.)  D

Contact Lens

CL Base Curve  D

CL Power  D

CL Over-Refracton (S.E.)  D

Note: Vertex Distance = 12.0

Outputs (D)

Calculate OS      Method

**OS**

Inputs

Pre Refractive Surgery

nc 1.3375

K1  K2

Spectacle Refraction (S.E.)  D

Post Refractive Surgery

nc 1.3375

K1  K2

Spectacle Refraction (S.E.)  D

Contact Lens

CL Base Curve  D

CL Power  D

CL Over-Refracton (S.E.)  D

Note: Vertex Distance = 12.0

Outputs (D)

Calculate      OK      Cancel

The best method in any instance will depend largely on the available patient information. The following section gives a brief description of the post-refractive formulas and their required inputs.

**History Derived** – The History Derived method is considered to be one of the most accurate methods for determining the true post-surgical corneal power. Unfortunately, it also requires the most information about the patient – some of which may not be available. This History Derived method calculates the amount of correction that was achieved through refractive surgery and adds it to the pre-surgical corneal power. This method can be used with 3<sup>rd</sup> generation formulas (SRK/T, Hoffer Q or Holladay 1).

#### Required Inputs

- Pre Refractive Surgery K measurements
- Pre Refractive Surgery Spectacle Refraction (S.E.)
- Post Refractive Surgery Spectacle Refraction (S.E.)

**Clinically Derived** - The advantage of this formula is that no information prior to the refractive surgery is needed. The only input required for this method is the average (post-surgical) manual keratometry measurement of corneal power. The disadvantage is that since the amount of correction achieved through refractive surgery is unknown, it is not factored into the calculation. This method can be used with 3<sup>rd</sup> generation formulas (SRK/T, Hoffer Q or Holladay 1).

#### Required Inputs

- Post Refractive Surgery K measurements

**Refraction Derived** - This method is useful when the pre-operative K values are unknown, but the pre-operative spectacle refraction is known.

#### Required Inputs

- Pre Refractive Surgery Spectacle Refraction (S.E.)
- Post Refractive Surgery K measurements
- Post Refractive Surgery Spectacle Refraction (S.E.)

**Contact Lens Over-Correction (CL Over-Corr)** - This method requires the doctor to place a plano hard contact lens of a known curvature on the patient's eye. If the refraction remains the same, the corneal curvature is equal to that of the contact lens. If the over-refraction results in a higher myopic error, the difference in refraction is added to the base curve value. This method can be used with 3<sup>rd</sup> generation formulas (SRK/T, Hoffer Q or Holladay 1).

Required Inputs

- Contact Lens Base Curve
- Contact Lens Power
- Contact Lens Over-Refraction (S.E.)
- Post Refractive Surgery Spectacle Refraction (S.E.)

**Double-K (SRK-T)** – This method uses the pre-refractive corneal power to estimate the post-refractive anterior chamber depth (also known as effective lens position). The post-refractive corneal power is estimated using the History Derived Method. This method can be used with the SRK/T formula.

Required Inputs

- Pre Refractive Surgery K measurements
- Pre Refractive Surgery Spectacle Refraction (S.E.)
- Post Refractive Surgery Spectacle Refraction (S.E.)

**\*\*\* IMPORTANT \*\*\***

The corrected corneal power (Kcorr) calculated using the formulas described above should *NOT be interpreted as a suggested or recommended corneal power on the part of DGH Technology, Inc.* The calculated corrected corneal power (Kcorr) should be used as a guideline only. It is based upon the accuracy of the post refractive formula used and the data entered. The determination of the corrected corneal power (Kcorr) to use for performing IOL Power Calculations must be made by the surgeon, based on personal experiences and previous postoperative results.

## 16.17 Using the Post Refractive Calculator



1. Select the **Post Refractive** checkbox on the IOL Calculator Screen next to the eye that has previously undergone corneal refractive surgery. This will automatically open the Post Refractive K Calculator. 

**Note:** The Post Refractive K Calculator can be opened again at any time by selecting the **Kcorr** button on the IOL Calculator Screen. The “Kcorr” button is inactive until the “Post Refractive” checkbox has been selected.




2. When the **Post Refractive** checkbox is selected, the Pre-Op group box on the IOL Calculator Screen will change in the following ways:
  - K1 and K2 will be replaced with Kcorr or Pre Ref. K and Post Ref. K (Pre Refractive Surgery and Post Refractive Surgery) depending on the Post Refractive Method selected.
  - The values Kcorr, Pre Ref. K, and Post Ref. K are “read only”
  - K1 and K2 are entered automatically into the “Post Refractive Surgery” section of the Post Refractive Calculator

**Note:** Deselecting the “Post Refractive” checkbox will return all fields to their original configuration.

3. In the Post Refractive K Calculator, select the desired eye(s) on which to perform Post Refractive calculations using the **Calculate OD** and **Calculate OS** checkboxes. 
4. Select the desired Post Refractive formula using the **Method** drop-down box. The required inputs for the selected method will then be highlighted in white. 
5. Enter all known inputs. If a field is gray, that input is not required for the selected method, but it can still be edited. K1 and K2 values can be entered as either mm or D, but both K values in a section must be in the same units. Pre-Surgical and Post-Surgical K values can be entered in different units.

**Note:** The refractive index (nc) displayed in the IOL Calculator will be automatically transferred to both the Pre-Surgical and Post-Surgical sections of the Post Refractive Calculator. If a different keratometer was used for Pre-Surgical measurements, the correct nc for that keratometer must be entered in the Post Refractive Calculator.

6. Select the **Calculate** button to calculate the corrected corneal power (Kcorr). If a required input is missing, the calculation will be cancelled and the necessary field will be highlighted in red. A successful calculation will display the Kcorr in the “Outputs” section of the calculator. Kcorr is given in diopters (D) using a keratometer index (nc) of 1.3375. 
7. Once calculated, the results will automatically be imported into the IOL Calculator Screen and used for IOL Power Calculations. Selecting the “OK” button without calculating will close the calculator and save the values entered, but will not import a Kcorr to the IOL Calculator.
8. Use the **K Correction Method** drop-down box on the IOL Calculator Screen to see the corrected corneal power (Kcorr) predicted by the other Post Refractive Formulas. If the new method selected requires information that has not been entered, the Post Refractive K Calculator will open automatically. 
9. Once the desired Post Refractive Formula has been selected, proceed with the IOL Power Calculation.

---

**Note:** Only IOL Power Calculation Formulas that are compatible with the currently selected Post Refractive Calculation Method will be available for selection in the IOL Calculator.

---

## 16.18 Performing IOL Calculations

The general procedure for performing IOL Calculations is as follows:

1. Select the **Patient Data** Screen and either enter a new patient or select an existing patient. Complete the fields in the Patient Info section of the page as described in Section 13.
2. Select the A-Scan Screen and perform either an Immersion or Contact A-Scan Measurement as described in Sections 14.9 and 14.10. A saved file from a previous measurement can also be used.
3. Select the IOL Calculator Screen. The measurement and patient data information will automatically be imported from the A-Scan and Patient Data Screens. Verify that the information is correct and enter any missing Pre-Operative Data.
4. Select the Doctor that is performing the Lens Replacement Surgery using the drop down box. The preferred IOLs and formulas will automatically be loaded for the selected Doctor.
5. The Scanmate software will automatically perform IOL Power Calculations once all necessary parameters for the subject eye(s) have been input.
6. If the Patient has previously undergone corneal refractive surgery, select the “Post Refractive” checkbox. This will automatically open the Post Refractive

Calculator. See Sections 16.16 and 16.17 for more information about performing Post Refractive Calculations.

7. Select the Measurement Source. Choices include:

- Unit Average
- Unit Single
- User Input

If User Input is selected, complete the required fields.

8. If desired, modify the IOLs and formulas using the appropriate drop-down boxes. The “default” lens constant for the selected IOL and formula will be shown underneath the formula selection box.
9. If desired, select the “Personalized Constant” checkbox. This option will only be available if a personalized lens constant has been calculated for the currently selected Doctor and IOL. When selected, the personalized lens constant will replace the default lens constant. (See Section 12.6 for more information on calculating personalized lens constants)
10. When all the inputs are correct, the IOL Calculator will automatically perform calculations and display the output for one or both eyes.
11. Save the pre-operative data for the current patient and procedure by selecting the **Save** button.

**\*\*\* IMPORTANT \*\*\***

The calculated IOL powers are *NOT to be interpreted as suggested or recommended implant powers on the part of DGH Technology, Inc.* The calculated IOL should be used as a guideline only. It is based upon the accuracy of the IOL formula used and the data entered. The final IOL implant power selection must be made by the surgeon, based on personal experiences and previous post-operative results.

## 16.19 Overview of IOL Power Calculation Formulas

The following IOL Power Calculation formulas are included in the DGH Scanmate IOL Calculator.

- SRK®-II
- SRK®-T
- Haigis
- Hoffer® Q
- Holladay 1
- Binkhorst II

The formulas selected for inclusion are commonly used, and have been refined from earlier models. IOL power calculation formulas have evolved from the first models to become more accurate in a wider range of cases; selecting the appropriate IOL formula will be based on the individual practitioner's judgment and experience.

There are two types of first generation formulas, theoretical and regression. The theoretical formulas, such as those developed by Thijssen, Colenbrander, Fyodorov, Binkhorst and van der Heijde are based on a two lens system, where the cornea and pseudophakic lens focus images on the retina. The first generation regression formulas were empirically derived from a large sample of patients that had undergone lens replacement surgery. The SRK® formula developed by Sanders, Retzlaff and Kraff is an example of a first generation regression formula. First generation IOL Power Calculation formulas calculated IOL power as a function of the measured Axial Length and Corneal Power. These formulas typically assumed a constant value for the post-operative anterior chamber depth. Due to this assumption, they often predict incorrect replacement lens powers for eyes that are longer or shorter than average.

The second generation formulas such as Binkhorst II, Hoffer®, Shammas and SRK®-II attempt to correct the errors in predicting replacement IOL powers that typically occur for long and short eyes. They did this by relating the predicted post-operative ACD to the axial length. The shortcoming of this approach is that it assumes that shorter eyes always have short anterior chamber depths and vice versa.

Third and fourth generation formulas have further refined the method for predicting the post-operative ACD. The SRK®/T, Hoffer® Q and Holladay 1 formulas predict the post-operative ACD using the measured corneal power (or radius) as well as the axial length. The Haigis formula uses the pre-operatively measured ACD and axial length to predict the post-operative ACD. All third and fourth generation formulas allow doctors to personalize their lens constants based on their surgical results with a particular IOL.

## **16.20 Factors Impacting the Accuracy of IOL Power Calculations**

Many factors play a role in predicting the correct IOL power for a patient. Several of the most common factors are described in the following sections.

### **16.20.1 Axial Length**

Axial length measurement is the most important factor affecting IOL power calculations. An error of 1.0 mm affects the postoperative refraction by approximately 2.5 diopters.

### **16.20.2 Corneal Power**

Corneal power is the second most important factor affecting IOL power calculations. A keratometer is an instrument that measures the central 3.3 mm of the anterior curvature of the cornea in its two meridians. The readings are called K-readings.

There are two potential sources of error in keratometry. First, failure to calibrate the instrument can cause all readings to be in error by as much as 0.2 mm or 1.0 D.

A second source of error is hidden in the diopter calculation of any keratometer. Keratometers cannot directly measure the refracting power of the cornea in diopters. Rather, they measure the radius of curvature of the anterior corneal surface and then convert this millimeter measurement into diopters. Because only the anterior surface of the cornea is measured, the effect of the posterior surface on the overall refracting power must be taken into account to convert the measurement into diopters. The dioptric scale used for this conversion is based on an adjusted refractive index. The true refractive index of the cornea is 1.376, but to estimate the true refracting power of the cornea from only the radius of the anterior surface, a fictitious refractive index ( $n_c$ ) is used. This index varies with the make of keratometer.

A radius of 7.8 mm will thus read 43.27, 43.08, or 42.56 D respectively, depending on which keratometer is used. These differences may not be of much importance in contact lens fitting or the determination of corneal astigmatism. However, these differences are significant when calculating IOL power.

Any IOL power calculation formula which requires the dioptric power of the cornea is subject to this source of error. With the same data, one may calculate IOL powers varying by almost a full diopter depending on the make of the keratometer used.

With both potential sources of error in mind, consider that an error of 1 diopter in the measurement of the corneal power produces an error of about 1 diopter in the postoperative refraction.

### **16.20.3 Postoperative Anterior Chamber Depth**

The anterior chamber depth is another factor that affects IOL power calculations. An error of 1.0 mm affects the postoperative refraction by approximately 1.0 diopter in a myopic eye, 1.5 diopters in an emmetropic eye, and up to 2.5 diopters in a hyperopic eye.

### **16.20.4 Surgical Technique**

Changes in corneal curvature are often noted postoperatively. This fact as well as differences between actual placement and predicted placement of the IOL can produce an error. An intraocular lens placed in the posterior segment requires a stronger power, and inserting the implant with the convex side backwards necessitates an even stronger lens.

### **16.20.5 Implant Power**

The implant power is measured differently by different manufacturers. This may require adjustments to IOL power predictions.

## 17 Creating Reports

Reports can be created for viewing from the Menu Bar. The DGH Scanmate software can produce IOL Calculator Reports, A-Scan Short Reports, A-Scan Custom Reports and B-Scan Reports. Reports can be printed out, saved as PDF files and added to the patient record database. All reports are created with a header that includes the Patient's Name, ID Number, Doctor and Operator. The revision number of the DGH Scanmate software used to generate the report is located in the footer of the report.

### 17.1 IOL Calculator Report


The IOL Calculator Report presents the results of calculations by the three selected IOL models for each eye. To create an IOL Calculator Report, use the menu bar at the top to go to **Reports → IOL Calculator Report**. A report will be created of the data currently available; any fields which are missing will appear as blank in the report. The standard IOL Calculations Report includes data for both OS and OD, and for three IOL models.

### 17.2 A-Scan Short Report

The A-Scan Short Report presents images of the waveforms captured and a table of measurements for each eye. To create an A-Scan Short Report, use the menu bar at the top to go to **Reports → A-Scan Short Report**. If there is no A-Scan data available, the software will prompt the user to create a scan.

### 17.3 A-Scan Custom Report

The A-Scan Custom Report allows the user to select which information to include in the report. To create a Customized Report, use the menu bar at the top to go to **Reports → A-Scan Custom Report**. A dialog box will open allowing the user to customize the report being created. Once the desired options have been selected, click on the "Generate Report" button to create the report.



The screenshot shows a dialog box for customizing an A-Scan report. It is organized into several sections with radio buttons and checkboxes:

- Eye:** Radio buttons for OD, OS, and Both. 'Both' is selected.
- Measurement Data:** Checkboxes for 'Measurement Bank' and 'Velocities', both of which are checked.
- Patient Fields:** Checkboxes for 'Additional Patient Info', 'Pre-Op Data', and 'Post-Op Data', all of which are checked.
- Waveform Data:** Radio buttons for 'All', 'Selected', and 'None'. 'All' is selected.
- IOL Calculations:** Radio buttons for 'All', 'IOL1', 'IOL2', 'IOL3', and 'None'. 'All' is selected.

At the bottom right of the dialog box are two buttons: 'Generate Report' and 'Cancel'.

## 17.4 Using Reports

All types of reports can be printed, browsed, and saved in the same ways.




Clicking the small printer icon in the top left corner will send the report to the printer.



The arrows and magnifier icon next to the printer icon will navigate through the report.



The **Export to PDF** button will save the report as a PDF to the Default Data Directory specified in System Preferences. A browser window will open to allow the user to select the save location and name of the PDF created. PDFs can also be exported as explained in Section 17.6.




The **Save to DB** button will add the report to the DGH-Scanmate database. Once saved, the report can be retrieved by searching the patient's name or ID number.



Clicking **Cancel** will exit the report and return to the Scanmate application.

## 17.5 Opening Reports

To view a previously saved report, select the patient through the Patient Data screen search. Using the Menu Bar, go to **File → Open** or click the **Open File** button.



A dialog box will display all the saved A-Scan measurements, videos and reports available for that patient. Each report is automatically given a unique name consisting of the type of report, the date it was created, and the time it was created.

## 17.6 Exporting A-Scan PDF Reports

It is possible to convert A-Scan reports saved in the Scanmate database to PDF files so that they can be shared with EMR systems or practitioners that do not have the Scanmate Software. To do this go to **File → Open** or click the **Open File** button.

This will open a window that shows all saved images, videos and reports for the current patient. Left-click and drag the report thumbnail(s) to the desired location. A PDF version of the report will be automatically generated.

## 18 Database Management

Patient data, scan images, measurements, video files, and reports created by the Scanmate application are saved in an DGH-Scanmate database. The database allows patient records to be shared, centralized, or accessed remotely, according to the needs of each biometry department. For example, several operators in a large practice could use several Scanmate Flex, DGH 6000 or DGH 8000 units simultaneously, but all patient records would be stored in a central database. Or, a doctor could perform scans in an examination room, then later access patient records from an office computer.

Refer to the Scanmate Installation Guide for information on backing-up and restoring, migrating, or deleting a database.

### 18.1 Importing and Exporting Data

Individual measurement and video files can be imported or exported from the DGH-Scanmate database to be consulted in other locations or by other practitioners. The Scanmate software must be installed on a computer to view measurement or video files.

To export a measurement file, use the menu bar at the top to go to **File → Export → A-Scan Measurements**. The measurements currently in the measurement bank will be saved as an .ames file.

To export a Scanmate video file, use the menu bar at the top to go to **File → Export → A-Scan Video**. The images currently in the video buffer will be saved as an .avdo file. Alternatively, an .avi file (which can be opened on most PC media applications) can be exported by using the top menu bar and selecting **File → Export → A-Scan Video (AVI)**.

To import a measurement file that is not currently in the database, use the menu bar at the top to go to **File → Import → A-Scan Measurements**. A browser window will open to select a .mes or .ames file to open. The opened file will replace any information currently in the patient data fields and the measurement bank. The opened file can be saved into the database at this point.

To import a video buffer that is not currently in the database, use the menu bar at the top to go to **File → Import → A-Scan Video**. A browser window will open to select a .vdo or .avdo file to open. The opened file will replace any information currently in the patient data fields and the video buffer. The opened file can be saved into the database at this point.

## 19 Electromagnetic Compatibility

Like other medical equipment, the DGH 6000 Scanmate-A requires special precautions to ensure electromagnetic compatibility with other electrical medical devices. To ensure electromagnetic compatibility (EMC), the DGH 6000 must be installed and operated according to the EMC information provided in this manual.

The DGH 6000 has been designed and tested to comply with EN 60601-1-2 requirements for EMC with other devices.



### CAUTION

Portable and mobile RF communications equipment may affect the normal function of the DGH 6000 Scanmate A.



### CAUTION

Do not use cables or accessories other than those provided with the DGH 6000 Scanmate A, as they may result in increased electromagnetic emissions or decrease immunity to such emissions.

### Guidance and Manufacturer's Declaration: Electromagnetic Emissions and Immunity

The DGH 6000 Scanmate A is intended for use in the electromagnetic environment specified below. The customer or the user of the DGH 6000 should ensure that it is used in such an environment.

Environmental Phenomena	Test In Accordance to	Level	Criteria	Basic Standard	Notes
Radiated Emissions	EN60601-1-2	Group 1 Class a	Under Limit	CISPR 11	Measure at 5 meters
Electrostatic Discharge	EN60601-1-2	±2Kv ±4Kv ±8Kv contact discharge  ±2Kv ±4Kv ±8Kv air discharge	36.202.1 (j)	EN61000-4-2	Apply to all accessible components
Radiated Immunity	EN60601-1-2	80MHz-2.5GHz 3V/m 80% @ 1kHz	36.202.1 (j)	EN61000-4-3	Expose all parts of EUT to field
EFT I/O Only	EN60601-1-2	±2Kv 5/50 5kHz	36.202.1 (j)	EN61000-4-4	None
Conducted Immunity I/O Only	EN60601-1-2	0.15 – 80MHz 3Vrms 80% @ 1kHz	36.202.1 (j)	EN61000-4-6	None

### **Guidance and Manufacturer's Declaration: Electromagnetic Immunity**

The DGH 6000 Scanmate A is intended for use in the electro-magnetic environment specified below. The customer or the user of the DGH 6000 should ensure that it is used in such an environment.

Field strength from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, armature radio, AM and FM radio broadcast and TV broadcast, cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the DGH 6000 is used exceeds the applicable RF compliance level, the DGH 6000 should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the DGH 6000.

## **20 Care and Maintenance**



### **WARNING**

Users of the DGH 6000 Scanmate-A have an obligation and responsibility to provide the highest degree of infection control possible to patients, co-workers and themselves. To avoid cross contamination, follow all infection control policies established for the office, department or hospital as they apply to personnel and equipment

### **20.1 Care of Transducer**

Be careful when handling the transducer. If the transducer is dropped on a hard surface it can be damaged. Regularly check the transducer front face for damage, as this may impair the performance of the device. The transducer may be submersed in water up to the cable (not wetting the connector) during normal use.

The transducer should be cleaned after every use, using the following protocol:

1. Wear protective gloves when performing the cleaning process.
2. Disconnect the transducer from the unit.
3. Use a soft cloth (or equivalent) dampened with a mild soap or 70% isopropyl alcohol solution to remove any particulate matter or body fluids that remain on the probe tip, probe housing and probe cable (if applicable). Do not immerse the entire probe in the cleaning solution. **Caution:** Care must be taken to prevent any liquids from coming into contact with the probe connector
4. To remove remaining particulates, rinse probe with water. Do not immerse the probe in the rinse solution.

5. Wipe with a dry cloth; or wipe with a water-dampened cloth to remove soap residue, and then wipe with a dry cloth.



**WARNING**

Do Not Attempt To Open the Transducer or USB Interface Module

## 20.2 Cleaning and Disinfecting the Transducer Tip

Keep the transducer tip clean and disinfected. To prevent patient-to-patient infection, a disinfection procedure should be followed after each patient. After disinfection, the tip should be rinsed in sterile distilled water before using.

The following disinfectants are recommended due to their compatibility with the probe tip material.

<b>Disinfectant</b>	<b>Concentration Tested*</b>
Cavicide Solution	(10-20%) Isopropyl Alcohol and (1-5%) Ethylene Glycol Monobutyl Ether
Cavicide Wipe	(10-20%) Isopropyl Alcohol and (1-5%) Ethylene Glycol Monobutyl Ether
Cidex	2.55% (w/w) Glutaraldehyde
Cidex OPA	6.2% by (w/w) Ortho-Phthalaldehyde (1,2 – benzenedicarboxaldehyde)
Isopropyl Alcohol	70% (v/v) Isopropyl Alcohol
Household Bleach	0.6% (w/w) Sodium Hypochlorite
Hydrogen Peroxide	3% (w/w) H <sub>2</sub> O <sub>2</sub>
Milton	2% (w/w) Sodium Hypochlorite

\*Note: The concentrations listed in this table are the specific concentrations that were tested by DGH to ensure compatibility with the probe tip material. DGH does not endorse or recommend the concentrations listed in the table above.

 **WARNING**

DGH makes no claims about the biological effectiveness as a disinfectant of any of the products listed above. Furthermore, DGH makes no claims regarding the effectiveness of any of these products for killing any known, or unknown, bacteria, virus, or other micro-organisms. DGH only claims that these products, when used properly, will not harm the transducer tip.

**Instructions for Disinfecting Transducer Surfaces when Using a Cavicide Wipe**

1. Wear protective gloves when performing the disinfecting procedure.
2. Detach the transducer from the unit.
3. Check the expiration date on the wipes that are being used. Use only wipes that are within the expiration date.
4. Thoroughly wet down the transducer tip, housing and cable (if applicable) according to the standard operating procedure (SOP) or equivalent document established by your facility or institution. **Caution:** Care must be taken to prevent any liquids from coming into contact with the probe connector.
5. Follow the instructions in the SOP for the duration of exposure to the disinfectant.
6. If applicable, follow the instructions outlined in the SOP for rinsing the probe tip.
7. Examine the probe for damage such as cracks, splitting, or sharp edges or projections. If damage is evident, discontinue use of the probe and contact a customer service representative.

**Instructions for Disinfecting Patient Contact Surfaces when Using a Liquid Disinfectant**

1. Wear protective gloves when performing the disinfecting procedure.
2. Detach the transducer from the unit.
3. Check the expiration date on the solution that is being used. Use only solutions that are within the expiration date.
4. Mix the disinfectant listed in the table above according to the standard operating procedure (SOP) or equivalent document established by your facility or institution.
5. Immerse the tip of the transducer approximately ¼ inch (6.35 mm) in the disinfectant solution.

6. Follow the instructions in the SOP for the duration of probe tip immersion. DGH probe tips may be damaged if immersed for longer than 1 hour.
7. If applicable, follow the instructions outlined in the SOP for rinsing the probe tip.
8. Examine the transducer for damage such as cracks, splitting, or sharp edges or projections. If damage is evident, discontinue use of the probe and contact a customer service representative.

 **WARNING**

It is the responsibility of the user to remain current with the latest information from the relevant disinfectant manufacturer concerning instructions, effects, necessary concentrations, immersion times and rinse requirements.

 **WARNING**

Using a non-recommended disinfectant, incorrect solution strength, or immersing the probe tip deeper than described in step #5 (above), or for a period longer than 1 hour can damage or discolor the probe tip and will void the probe warranty.

Do not immerse the probe tip for longer than one hour. The probe tip may be damaged by longer immersion times.

Disinfect the probe tip using only the liquid solutions above. Using autoclave, gas (EtO), or other non DGH Technology approved methods will damage the probe and void the warranty.

 **WARNING**

The transducer should **NEVER** be autoclaved or subjected to intense heat. As a general rule, the above cleaning instructions are sufficient to disinfect the transducer in ordinary use. Do not scratch or chip the conical transducer tip, which makes contact with the cornea.

### **20.3 Cleaning and Disinfecting the Immersion Shell**

The immersion shell is not sold sterile, and must be cleaned and disinfected before each use. The immersion shell cannot be autoclaved. Clean as described in the immersion shell Instructions for Use (included in package). The catheter tubing cannot be reused.

### **20.4 Care of USB Interface Module**

Keep the USB transducer plug dry at all times. Regularly check the transducer and USB cables for cuts, cracks and kinks. The presence of these defects can impair the performance of the device.

### **20.5 Cleaning the USB Interface Module**

The DGH 6000 USB Interface module housing and USB cable can be cleaned using a cloth dampened with a mild soap and water solution. The connectors should not be exposed to water.

### **20.6 Operating Conditions**

The DGH6000 (Scanmate A) should be operated between temperatures of 18°C and 40°C.

### **20.7 Verifying A-Scan Calibration**

The DGH 6000 uses a pattern recognition algorithm that looks for unique characteristics of the human retina to obtain true probe alignment. Under normal usage, the DGH 6000 will not recognize standard plastic interface phantoms as eyes to be measured. For this reason, the DGH 6000 must be placed into a test mode to check calibration on a standardized polystyrene block, provided with the unit. The procedure for entering the test mode and verifying proper calibration of the unit is described below.

1. Navigate to the Patient Data Page by selecting the “Patient Data” tab. In the “Last Name” drop-down box, select “Test Block”.
2. Navigate to the A-Scan Page by selecting the “A-Scan” tab. The current eye should be displayed as “Test Block” in the waveform display screen, and the lens type drop-down box should be frozen, displaying “Test Block”.
3. Set the provided polystyrene reference block on a stable, level surface. Apply a small drop of water to the end of the reference block.
4. Applanate the probe tip to the center of the reference block. The probe tip and reference block have equal diameters to facilitate positioning.

5. The DGH 6000 should measure the reference block almost instantaneously. If not, move the probe tip slightly until a measurement is obtained.
6. Observe axial length measurements to be within the stated tolerance marked on the reference test block.
7. Readings outside these limits do not necessarily indicate the DGH 6000 is out of calibration. Erroneous measurements may be due to the following:
  - a) During the test mode, the algorithm does not preclude measurement of tear film as it does during normal operation. Therefore, if excessive water was placed on the reference block, the water depth may be measured causing longer axial length measurements.
  - b) The polystyrene reference block is sensitive to temperature. The block must be stabilized at  $72^{\circ}\text{F} \pm 10^{\circ}\text{F}$  before measurements are taken.

If an erroneous reading is obtained, verify that these fault conditions are not present and repeat the measurements on the reference block.
8. If the reference block readings are outside the limits stated above, it is possible that the unit is defective and needs repair. Contact the Customer Service Department of DGH Technology, Inc. at (610) 594-9100 to arrange for repair.
9. If desired, the results of the calibration check can be saved in the “Test Block” patient folder.

## 20.8 Storage

When the DGH 6000 is not being used, it should be stored in a clean, dry area.

To prevent damage to the DGH 6000, do not store in areas where it might be exposed to:

- Excessive vibration
- Excessive dust and dirt
- Liquids or condensation
- Impact

Store the DGH 6000 under the following ambient conditions:

- Temperature:  $-10^{\circ}\text{C}$  to  $50^{\circ}\text{C}$  ( $14^{\circ}\text{F}$  to  $122^{\circ}\text{F}$ )
- Relative Humidity: 20% to 80%
- Atmospheric pressure: 70 kPa to 106 kPa

## **20.9 Transportation**

Never carry the DGH 6000 by the USB or transducer cables.

Never bend the USB or transducer cables in a tight radius. This could result in damage to the cables.

Transport the DGH 6000 under the following ambient conditions:

- Temperature: -10°C to 50°C (14°F to 122°F)
- Relative Humidity: 20% to 80% (no condensation)
- Atmospheric pressure: 70 kPa to 106 kPa

When transporting the DGH 6000 to a different field location or when returning it for repair or maintenance, use the original DGH 6000 packing enclosure.

If the original package is not available, pack in such a way that the DGH 6000 is protected.

## **20.10 Disposal**

Contact DGH Technology, Inc. before disposing of the DGH 6000.

Concerning the WEEE label, the following information is for EU member states:

The use of this symbol indicates that this product should not be treated as household waste. By ensuring that this product is disposed of correctly, you will help prevent potential negative consequences for the environment and human health, which could otherwise be caused by inappropriate waste-handling of this product. For more information concerning the return and recycling of this product, please consult DGH Technology, Inc.

## **21 Troubleshooting**

Refer to the Scanmate Installation Guide for troubleshooting instructions related to the installation and configuration of the Scanmate software. The installation guide also contains information on backing up, restoring and moving the patient database.

### **21.1 Unable to Get a Measurement**

- Verify that the correct measurement mode (Immersion/Contact) and “Eye Type” has been selected.
- If performing a Contact measurement, it may be necessary to adjust the Compression Sensitivity (Comp). See Section 14.8 for details.
- If performing an Immersion measurement, verify that the transducer is properly inserted in the Prager Shell. See Section 14.10 for details.
- There will be some eyes which cannot be measured in Automatic Measurement (Auto Meas) mode. In these cases, measurements can be made by using the Manual Measurement (Manual Meas) mode.

## **22 Warranty**

DGH Technology, Inc. “DGH” warrants each new DGH 6000 and its accompanying accessories (hereinafter called “Equipment”) to be free from defects in material and workmanship for twelve (12) months from the date of delivery to the original purchaser. This warranty is not applicable to any defect that is the result of an accident, misuse, mishandling, neglect, improper installation, improper repair or improper modification by persons other than DGH. This warranty does not apply if the Equipment has not been operated and maintained in accordance with the operating and maintenance manuals and instructions or bulletins issued in respect thereof by DGH. The cost of servicing replaceable and expendable items including parts and labor made in connection with the routine maintenance services as described in such Operator’s Manual is not covered under this warranty and is the responsibility of the purchaser.

This warranty is strictly limited to replacement or repair of the part that is found to be defective in material and workmanship. At the option of DGH, said part shall be replaced or repaired free of charge, F.O.B. our factory by DGH.

DGH reserves the right to make changes in the design and material of Equipment without incurring any obligations to incorporate such changes in Equipment already completed on the effective date of any such changes.

This is the only warranty of this product and is expressly in lieu of all other warranties, expressed or implied by law or otherwise, including any implied warranties of merchantability and of fitness for a particular purpose. Without regard to the alleged

defect, DGH does not, under any circumstances, assume any responsibility for the loss of time, inconvenience or other consequential damages, including but not limited to, loss or damage of personal property, or loss of revenue. DGH has neither assumed nor authorized any other person (including any distributor authorized to sell its Equipment) to assume for it any other liability in the connection with the sale of Equipment.

## 23 Lifetime / Shelf-life

The shelf-life / lifetime indicated for this device is 10 years.

## 24 Customer Service

If you are having problems with this unit, please refer to the appropriate sections of this manual. Most service calls result from a misinterpretation of the operation of the instrument. The instructions in this manual have been carefully reviewed to ensure error-free performance of the DGH-6000.

However, if you feel there is a problem with the unit or a transducer, please contact the Customer Service Department at the address below. DGH Technology, Inc. can also be contacted via our website at [www.dghtechnology.com](http://www.dghtechnology.com). When contacting us, please provide the model and serial number for the unit. The model number and serial number are located on the underside of the USB Interface module and can also be viewed on the display by selecting the “About” button found on the “Help” toolbar.

Manufactured by:

**DGH TECHNOLOGY, INC.**  
 110 SUMMIT DRIVE  
SUITE B  
EXTON, PA 19341  
USA (610) 594-9100

**DGH** *TECHNOLOGY, INC.*  


Authorized European Representative:

### **EMERGO EUROPE**



Prinsessegracht 20  
2514 AP, The Hague  
The Netherlands

## APPENDIX A COMPUTER SYSTEM SPECIFICATIONS

### Minimum System Requirements

<b>Processor:</b>	32 bit or 64 bit, 2 GHz
<b>Memory:</b>	2 GB RAM
<b>Hard Drive:</b>	1 GB minimum, 100 GB recommended
<b>Ports:</b>	USB 2.0
<b>Display:</b>	1024 x 768 Resolution
<b>Peripherals:</b>	Mouse (or Touchpad), Keyboard and Speakers (For Alignment Audible Feedback)
<b>AC Power Supply:</b>	Medical Grade

### Compatible Operating Systems

Microsoft Windows 7 (32-bit or 64-bit)	Microsoft Windows Server 2008 R2 (64-bit)
Microsoft Windows 8/8.1 (32-bit or 64-bit)	Microsoft Windows Server 2012 /2012 R2 (64-bit)
Microsoft Windows 10 (32-bit or 64-bit)	Microsoft Windows Server 2016 (64-bit)

### WARNING

The use of a “Non-Medical” grade AC Adapter could potentially cause harm to the system, the operator and/or the patient.

### WARNING

Using “Non-Essential” Software in Conjunction with the Scanmate System Could have Unknown / Adverse Impact on the Operation of the Device and is Therefore Not Recommended.

### WARNING

Due to the Threat of Computer Viruses, it is Recommended that an Anti-Virus Program be Installed on the Computer Running the Scanmate Application and that Patient Records Be Backed up Regularly.

**APPENDIX B      SCANMATE-A SPECIFICATIONS**

Transducer	Starting Frequency: 12.5 MHz +/- 1.25 MHz. Damped Frequency: 10.0 MHz (nominal) Fixation LED Focused Acoustic Beam (23.0 mm nominal) Detachable (via Lemo Connector) Probe Dimensions: L 55.7 mm x D 6.6 mm Cord Length: 195 cm
USB Interface Module	Dimensions: L 146 mm x W 88 mm x H 38 mm Weight: 285 gm Interfaces: USB 2.0 + Lemo 4 pin Connector Power Requirements: 5 VDC, 500 mA (2.5 W max). Integrated Probe Holder
Measurement Range	Axial Length (AXL): 15 mm to 40 mm Anterior Chamber Depth (ACD): 2.0 mm to 6.0 mm Lens Thickness (LT): 2.0 mm to 7.5 mm
Accuracy	Measurement Repeatability: $\pm 0.03$ mm STDEV (Immersion) Resolution: 0.01mm
Measurement Modes	Water Immersion Contact
IOL Formulas	SRK-T, SRK-II, Binkhorst II, Holladay 1, Hoffer Q, Haigis
Post Refractive Formulas	Double K, History Derived, Clinically Derived, Refraction Derived, Contact Lens
Archive Functions	Patient Data, Exam Information, Measurement Files, Video Files, IOL Calculation Reports
Environmental	Operating Temp 18°C to 40°C Storage Temp -10°C to 50°C Relative Humidity 20% to 80% (no condensation) Atmospheric pressure 70 kPa to 106 kPa

## APPENDIX C      REFERENCES

IOL formulas: The IOL Power Calculations are based on formulas published in peer-reviewed journals.

### Binkhorst II

JEDMED Internal documents, based on Binkhorst, R.D. "Pitfalls in Determination of Intraocular Lens Power Without Ultrasound." *Ophthalmic Surgery* 7:68 (1976): 69-82.

STORZ CompuScan LT. User Manual.

### Haigis

Haigis, W. Personal communication, based on documents prepared July 2, 2008. Received Feb 2, 2009.

### Hoffer Q

Hoffer, K.J. "The Hoffer Q formula: A comparison of theoretic and regression formulas." *Journal of Cataract and Refractive Surgery* 19:6 (Nov 1993): 700-712.

Errata in *Journal of Cataract and Refractive Surgery* 33:1 (Jan 2007): 2.

Author response in *Journal of Cataract and Refractive Surgery* 33:1 (Jan 2007): 2-3.

### Holladay 1

Holladay J.T., et al. "A three-part system for refining intraocular lens power calculations." *Journal of Cataract Refractive Surgery* 14:1 (Jan 1988): 17-24.

### SRK-II

Sanders, D.R., et al. "Comparison of the SRK-II Formula and Other Second Generation Formulas." *Journal of Cataract and Refractive Surgery* 14:2 (Mar 1988): 136-141.

### SRK-T

Retzlaff, JA et al. "Development of the SRK/T Intraocular Lens Implant Power Calculation Formula." *Journal of Cataract and Refractive Surgery* 16:3 (May 1990): 333-340.

Errata in *Journal of Cataract and Refractive Surgery* 16:4 (Jul 1990): 528.

## Post Refractive Surgery Formulas

By in large, these methods have been implemented as described in H. John Shammas's book "Intraocular Lens Power Calculations" published by Slack Incorporated 2004.

### History Derived

$$K_{\text{corr}} = \text{pre-}K_{\text{avg}} - \left[ \frac{\text{S.E.}_{\text{post}}}{(1 - \text{S.E.}_{\text{post}})} - \frac{\text{S.E.}_{\text{pre}}}{(1 - \text{S.E.}_{\text{pre}})} \right]$$

### Refraction Derived

$$K_{\text{corr}} = \text{post-}K_{\text{avg}} - 0.23 \times \left[ \frac{\text{S.E.}_{\text{post}}}{(1 - \text{S.E.}_{\text{post}})} - \frac{\text{S.E.}_{\text{pre}}}{(1 - \text{S.E.}_{\text{pre}})} \right]$$

### Clinically Derived

$$K_{\text{corr}} = (1.14 \times \text{post-}K_{\text{avg}}) - 6.8\text{m}^{-1}$$

### Contact Lens Over Correction

$$K_{\text{corr}} = \text{CL}_{\text{base curve}} + \left[ \frac{\text{CL}_{\text{over}}}{1 - 0.012\text{m} \cdot \text{CL}_{\text{over}}} \right] - \left[ \frac{\text{S.E.}_{\text{post}}}{(1 - \text{S.E.}_{\text{post}})} \right] + \text{CL}_{\text{power}}$$

### Double K

The Double K method modifies the SRK-T formula to use both pre- and post- refractive surgery K values. The pre-refractive surgery K value is used to calculate the effective lens position, and the post-refractive surgery K value is used for IOL power calculations.

Aramberri, J. "Intraocular Len Power Calculation After Corneal Refractive Surgery: Double-K Method." *Journal of Cataract and Refractive Surgery* 29:11 (Nov 2003): 2063-2068.

## APPENDIX D TERMS AND ABBREVIATIONS

The following terms and abbreviations are used within this Operator's Manual and in the Scanmate program:

A-constant	The A-constant is a constant provided by an IOL manufacturer to be used in IOL power calculations.
ACD	Anterior Chamber Depth is defined as the measured distance from anterior vertex of the cornea and the anterior vertex of the lens.
ACD (constant)	The ACD constant is a constant provided by the IOL manufacturer to be used in Binkhorst formula IOL power calculations.
a0	The a0 constant is calculated from the manufacturer's A-constant to be used in Haigis formula IOL power calculations.
AXL	Axial Length is the distance between the front face of the cornea and the front face of the retina of the eye.
CL <sub>base curve</sub>	The contact lens base curve is the known curvature of the plano hard contact lens.
CL <sub>over</sub>	The contact lens over correction is the measured refraction with the plano lens in place.
CL <sub>power</sub>	The contact lens power is the known power of the plano hard contact lens.
IOL	Intra Ocular Lens
K <sub>corr</sub>	The corrected K value that should be used in IOL calculations for patients who have previously undergone refractive surgery.
K1	Corneal Power in Vertical Meridian (D)
K2	Corneal Power in Horizontal Meridian (D)
LT	Lens Thickness
nc	The keratometer constant.
pACD	The pACD constant is calculated from the manufacturer's A-constant to be used in Hoffer Q formula IOL Power calculations.
post-K <sub>avg</sub>	Average of K1 and K2 values measured after refractive surgery.
pre-K <sub>avg</sub>	Average of K1 and K2 values measured before refractive surgery.
Rx	Refractive error (D)
S.E. <sub>post</sub>	Refractive error measured after refractive surgery.
S.E. <sub>pre</sub>	Refractive error measured before refractive surgery.
SF	The Surgeon Factor constant is calculated from the manufacturer's A-constant to be used in Holladay 1 formula IOL Power Calculations