

## HANDLE OPTIONS



**2 STEP: 2 SETTINGS + FULL EXTENSION**  
Ultra light titanium for LRI, tunnel, side port, clear cornea and depth controlled grooves. Accepts 1mm to 3.5mm blades.



**MICROMETER: DIAL IN DEPTH**  
Ultra light titanium for AK, radial and transverse incisions. Accepts all 1mm blades.



**FREESTYLE: 1MM-3.5MM STRAIGHT OR ANGLED**  
Ultra light titanium for clear cornea, tunnel, cataract, paracentesis. Accepts 1mm to 3.5mm wide blades.



**LRI: DUAL PLATE, PRESET DEPTH**  
Ultra light titanium for LRI. Footplates and blade "pop" into position as handle is engaged. Accepts 1mm blades.



**EXCEL SLIM**  
Ultra light titanium. Economical. Freestyle cuts or preset for groove, paracentesis.

Each **DGH-KOI™** knife is manufactured to precise specifications. Our 28 step quality assurance check guarantees excellence. All the **KOI™** handles are manufactured from titanium and have a lifetime warranty. **DGH-KOI™** supplies a full line of knives for your surgical procedures. For more information, call 800-344-9674 to speak directly with a technician regarding your specific needs.

**DGH-KOI™** specializes in the **manufacturing and repair** of surgical diamond knives.

Our **design team** is available to personally assist you in the selection of your diamond knife.

Inquire about our **exchange program**—the **DGH-KOI™** technicians are experts in the modification and updating of your existing instrumentation. The older blade styles you have invested in are a valuable resource and can be reconfigured to meet your current surgical needs. Take advantage of our **in-house blade honing and reconfiguration service**.

At **DGH-KOI™**, we take pride in providing the most comprehensive warranty in the industry.

**KOI™** began designing and building diamond knives in 1980. We have built our reputation on delivering the finest quality products at affordable prices.

As a pioneer and innovator, **KOI™** remains committed to providing the ophthalmic community with tomorrow's technology today.



**BOX 120  
5055 SPRING RD  
SHERMANSDALE, PA 17090**

**PHONE: 800-344-9674  
INTERNATIONAL: 717-582-3020  
FAX: 717-582-2406  
EMAIL: [INFO@DGHKOI.COM](mailto:INFO@DGHKOI.COM)  
[WWW.DGHKOI.COM](http://WWW.DGHKOI.COM)**



## DIAMOND KNIFE QUICK REFERENCE GUIDE



### **DGH-KOI™ The Diamond Knife Specialist**

Select your handle and blade configuration from the enclosed options.

If your particular blade configuration is not depicted in the design options, inquire about our specialized in-house blade design program. Our technical team will meet your specific needs.

**800-344-9674**



## FOOTPLATE


### OPTIONS:




LRI UNIVERSAL



SINGLE PLATE



UNIVERSAL DUAL



VISI PLATE

## FREESTYLE TIP

### OPTIONS:



1MM STRAIGHT



1MM-3.5MM ANGLE



2MM-3.5MM STRAIGHT

## BLADE CONFIGURATIONS

### TAPERED SPEARS



Standard taper:	
WIDTH	BLADE#
1.5/2.7MM	1658
2.1/2.4MM	1618
2.3/2.9MM	1668
2.5/3.0MM	1617
2.5/3.2MM	1671
2.5/3.5MM	1640
2.6/3.0MM	1641
2.7/3.0MM	1622
2.7/3.2MM	1609
2.7/3.5MM	1623
2.8/3.1MM	1634
2.8/3.2MM	1656
3.0/3.2MM	1636
3.0/3.5MM	1648
3.2/3.5MM	1630
3.3/3.5MM	1611
Stealth/Large bevel-up:	
2.5/3.0MM	1624
2.5/3.5MM	1639
2.7/3.2MM	1628

### SPEARS



Standard:	
WIDTH	BLADE#
1.0MM	1600
2.0MM	1602
2.5MM	1612
2.8MM	1603
3.0MM	1605
Sharp sides:	
1.0MM	1601
2.5MM	1616
2.65MM	1626
2.75MM	1641
2.8MM	1604
3.0MM	1606
3.2MM	1608
3.5MM	1615
Stealth/Large bevel-up:	
2.5MM	1621
2.65MM	1627
2.8MM	1625
3.0MM	1620
3.2MM	1632



Inquire about YOUR custom width!

### LANCES



Standard:	
WIDTH	BLADE#
1.0MM	1500
1.5MM	1505
2.0MM	1502
2.8MM	1508
3.0MM	1503
3.2MM	1504
Sharp sides:	
1.0MM	1511
2.0MM	1517
2.8MM	1507
3.0MM	1513
3.2MM	1519

### CRESCENTS



WIDTH	BLADE#
1.7MM	1004
2.0MM	1000

## DIAMOND KNIFE REPAIRS

DGH-KOI™ manufactures and repairs all of our diamond knives in-house. We provide 3 day turnaround on the repair of KOI™ knives and many competitor knives as well. The evaluation and quote are free.

800-344-9674

## SINGLE EDGES



WIDTH	BLADE#
1.0MM	1100



WIDTH	BLADE#
1.0MM	1102



WIDTH	BLADE#
1.0MM	1103

## DOUBLE EDGES



WIDTH	BLADE#
1.0MM	1200



WIDTH	BLADE#
1.0MM	1201



WIDTH	BLADE#
1.0MM	1203

## TRIPLE EDGES



WIDTH	BLADE#
1.0MM	1402



WIDTH	BLADE#
1.0MM	1404



WIDTH	BLADE#
.5MM	1302
.75MM	1305
1.0MM	1301



WIDTH	BLADE#
1.0MM	1300



WIDTH	BLADE#
1.0MM	1304



WIDTH	BLADE#
1.0MM	1405

## EXCEL SERIES



#1701  
SINGLE  
EDGE  
45°



#1703  
LANCE  
30°



#1700  
SINGLE  
EDGE  
30°



#1702  
TRIFACET  
30°

PHONE: 800-344-9674  
INTERNATIONAL: 717-582-3020  
FAX: 717-582-2406  
EMAIL: INFO@DGHKOI.COM  
WWW.DGHKOI.COM

